# F·S·T HANDCRAFTED PRECISION For your




FINE SCIENCE TOOLS
2026

### **CHECK OUT OUR NEW WEBSHOP**

#### finescience.com









#### **LETTER FROM THE OWNER**

Markus Gawenda

#### Dear Valued Customers, Partners, and Friends,

As we celebrate our 51st year of precision, passion, and purpose, we extend our deepest gratitude to you; our valued partners, clients, and colleagues who all are an integral part of this remarkable journey. From our humble beginnings to our global presence today, one thing has never changed: our unwavering commitment to quality, service, support, and innovation.

Behind every catalog page is a dedicated team that truly cares – not just about instruments, but about the people who use them. Whether you're a long-time customer or opening our catalog for the first time, we are here to guide you, collaborate with you, and deliver an experience rooted in trust, quality, and care.

This year, we are proud to unveil our theme: Smallville – a tribute to our unwavering attention to detail and profound commitment to craftsmanship. It serves as a powerful and imaginative symbol of the qualities that define us. Just as each miniature figure is meticulously crafting every cutting

edge, every curve, every tip and mechanism reflects the heart and soul we pour into our instruments. Every element matters and every customer deserves nothing less than excellence.

This commitment to precision extends beyond our instruments, it is also embedded in every aspect of our processes, on how best to serve you. We at Fine Science Tools (FST®) are always looking for smarter, more sustainable ways to help our community. We are actively rethinking materials, refining production processes, and protecting the health of our planet alongside the health of everyone. We are making conscious efforts at FST and every decision we make today considers the needs of tomorrow. Thank you for being a part of our story. Here's to the next chapter – built on precision, shaped by innovation, and guided by values that never go out of style.

Warmest regards, Markus Gawenda





#### PICK YOUR PERFECT match

We value the smallest details - including your preferences. Choose whether you'd like to receive the FST catalog in digital, physical, or both formats.

Please refer to the enclosed bookmark for all the steps, or scan the code here:




https://fst-gr.com/mycatalog

### **TIMELINE**

#### **FST® SINCE 1974**

For more than 50 years, Fine Science Tools has been the leading distributor of precision surgical and microsurgical instruments, as well as laboratory accessories with offices in Foster City (USA), Heidelberg (Germany) and Hong Kong (China) supporting and serving customers across the globe.



2024

### **CREATIVE DESIGNS**

#### **CATALOG COVERS AND PRINT ADS**


The unveiling of our annual marketing campaign reflects our passion and excitement for innovation and creativity year after year. FST® catalog covers and print Ads have won over 22 awards and are well known within the research community.















































































### **EXCLUSIVE LINE**

#### **BONN-GAWENDA SCISSORS**

Our Exclusive Line was released on our 50th Anniversary, these instruments are meticulously crafted to support every aspect of your research needs. At the heart of this collection is the Bonn-Gawenda Scissors, created as a tribute to our founder, Hans Gawenda, and symbolizing our legacy of innovation and excellence. These scissors stand out for their advanced design and exceptional craftsmanship, reflecting the high standards that have always defined our work.

The ToughCut® technology incorporates one micro-serrated cutting edge paired with a razorsharp blade, delivering consistently clean and precise cuts. These scissors have perfect tips for precision work under a microscope. Designed with comfort and ergonomics in mind, the large finger loops and specially engineered handle ensure exceptional control throughout each use.








Our Exclusive Line offers a unique combination of popular features. Whether you're seeking elegance, stability and ergonomics, or versatile functionality, the Exclusive Line was specially designed to cater to every aspect of your needs.



**Look for the Exclusive Line** icon in the catalog.



#### **TABLE OF CONTENTS** | CATALOG 2026

Letter from the Owner	1
Timeline	2 - 3
Creative Designs	4 – 5
Exclusive Line	6 - 7
Table of Contents	8 - 9
Catalog Features	10 – 11
Precision that makes a mark	12 – 13
Handcrafted, Custom Design	12
Quality Control, Scissors Test	13
Our Range of Services	14 – 15
Overview	14
Sharpening & Repair	15
Alloys & Materials	16 – 17
Maintenance & Care	18 – 21











	23 - 55
Spring	23 - 34
Fine	35 - 48
Surgical	49 - 55
	56 - 69
Rongeurs	56 - 59
Cutters	60 - 66
Accessories	67 - 68
Curettes & Chisels	69
	70 – 81
	Fine Surgical  Rongeurs Cutters Accessories

Forceps/Tweezers		83 - 115	
F	ST® by Dumont	83 - 92	
	Fine	93 - 96, 99 - 103	
	Moria	97 - 98	
	S&T®	98	
	Standard	104 - 115	
Hemostats		116 - 121	
Vascular Instruments		122 – 131	
	Clamps	122 - 126	
Occluders & Ligation Dilation & Cannulation		127 - 128	
		129 - 131	

Probes & Hooks	133 – 138
Spatulae	139 - 142
Spoons	143
Pins & Holders	144 – 145
Wound Closure	146 - 156
Needle Holders & Suture	146 - 153
Staplers, Clips & Applicators	154 - 156
Retractors	157 – 163

Animal Accessories	165 - 170
Feeding Needles	165 - 166
Animal Identification	167 – 170
Instrument Care & Accessories	171 – 182
Instrument Care	171 – 172
Sterilization & Storage	173 – 182
Surgical Accessories	183 – 186
Drill Accessories	183
Cauterizers	184 - 186
Laboratory Accessories	187 – 197
Tools	187 – 191
Magnifiers & Lens Paper	192 - 197

Student Quality Instruments	200 - 217
Quality instruments with the student in mind	
Scissors	200 - 205
Scalpel Handles	206
Forceps	207 - 212
Hemostats	213 - 215
Pin & Needle Holders	216
Retractors, Vessel Dilator & Clamps	217




Index Pages		218 - 225
	Numerical Index	218 – 221
4-7-1	Alphabetical Index	222 - 225

### **CATALOG FEATURES**

#### **PRODUCT FEATURES**

To provide you with our support regarding your decisions in finding the right instruments for your purpose, the product feature details are explained below.



#### **ICONS**



The Exclusive Line offers a unique combination of popular features of our tools combined into one. Whether you're seeking elegance, stability and ergonomics, or versatile functionality, the Exclusive Line is specially designed to cater to every facet of your needs.



The tips have micro serrations to ensure a better grip on tissue or other materials.



These items are not suitable to autoclave, but can be sterilized by a different method. For example: hot bead sterilization, gas sterilization, plasma sterilization or hot air sterilization. Please contact us to receive more information regarding the best preferred method for your instruments.



Fine Science Tools has noted which items are most frequently purchased for use in a variety of surgical procedures. To find our most popular items, look for the "most popular" icon throughout the catalog.



Our yearly selection of new instruments is based on a careful study and knowledge of our customer needs, as well as listening to their requests.



Fine Science Tools is proud to offer a large selection of left-handed surgical instruments to our community. To find our left-handed instruments in this catalog, look for the "lefthanded" icon on the pages below.

PRICES SUBJECT TO CHANGE PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER

Scissors: pages 38, 49 Hemostatic Forceps: page 117 Needle Holders: pages 150 - 151 Ligation Aid: page 128



FST® features a selection of scissors with larger than normal finger loops. Most of these loops are easily adjustable to fit larger fingers more comfortably. To find our large loop scissors please look for the "large loop" icon on the pages below.

Scissors: pages 35, 39, 48 and 51



#### **WE PROUDLY STOCK**

Fine Science Tools offers a wide selection of the latest and most popular Dumont forceps. Handcrafted in Switzerland for more than a century, these forceps are known for their consistent and uncompromised fineness.

Dumont forceps come in an array of styles and alloys. Below you will find detailed descriptions to help guide you in choosing the right Dumont forceps best suited for your needs.

#### FINDING THE DUMONT FORCEPS BEST SUITED FOR YOUR NEEDS

Alloys	Hardness	Magnetic	Temperature resistance	Autoclavable	Advantages
Titanium	Least Hard	No	Up to 430 °C	Yes	Resistant to corrosion from nitric acid, chloride, and salt water. Titanium is 40% lighter than Surgical Stainless Steel.
Dumoxel	Hard	No	400 °C	Yes	The most popular alloy available for Dumont forceps offers the high resistance to corrosion, sulphuric environments, hydrochloric acids, and all other mineral and organic acids.
lnox	Harder	Yes	400 °C	Yes	Available for most Dumont forceps. Contains chromium added to carbon steel, which results in some loss of hardness but offers good stainless qualities and resistance to corrosion.
Dumostar	Harder	No	500°C	Yes	More durable and corrosion resistant than the best Surgical Stainless Steel. Highly resistant to mineral and organic acids, and salt water. Dumostar lasts an average of four times longer than Dumoxel and Inox. Most cost effective alloy for laboratory use.
Carbon	Hardest	Yes	N/A	No	Extremely hard; ensuring strong tips. Easily stains and rusts without proper care. Carbon cannot be immersed in a water bath or sterilized.

#### Tip profiles: Dimensions vary. See each style for specific dimensions (width x thickness).

NOT FOR HUMAN USE PRODUCTS ARE NOT INTENDED FOR USE ON OR IN HUMANS - MOST ITEMS SHOWN IN ACTUAL SIZE


Standard	Made for high precision and consistently precise work under a microscope.
Biology	Twice as fine as "Standard Tips". Especially produced for high precision laboratory work under the microscope.
Super Fine Tips	Four times finer than the "Standard Tips". These tips are extremely delicate for the finest work under the microscope.
	Extra care is needed for these tips as they are very fragile.

#### Styles

,	
#2	Straight forceps with wide yet fine tips, referred to as Mouse Laminectomy Forceps because of their popularity in breaking off delicate bone from mouse skulls.
#2AP	Epoxy coated with flat wide blunt tips.
#3	Long straight forceps with coarse tips. Popular in low magnification procedures.
#3C	Shorter version of #3.
#4	Similar to the #5 but with broader shanks and stronger tips.
#5	Most popular style of Dumont straight forceps with fine tips which come in several styles.
#55	Similar to the #5 but with lighter shanks, giving the forceps a lighter spring action.
#5CO	Similar to the #55 but with even finer shanks producing the lightest spring action and longest meeting surface of any Dumont style.
#5L	Medical #5 forceps with slide sleeve for a secure hold and retention.
#5SF	The finest tips ever produced by Dumont.

#5XL	Extra Long #5 forceps.
#5/45	Tips angled at 45°.
#5/45C	Designed to pick up cover slips from multi-well plates.
#5/90	Tips angled at 90°.
#5AC	Self-closing Anti-Capillary forceps.
#6	Sharply angled, broad tips.
#7	Most popular style of Dumont curved forceps with fine tips which come in several styles.
#7B	Curved serrated tips.
AA	Strong broad tips.
М	Mini forceps.
Medical	Denotes distinctive ribbed handle to maximize grip.
SS	Long and narrow shanks with thin tips.
WA	One piece design with no soldering grooves allows the forceps to be thoroughly and easily cleaned preventing biofilm residue.

### PRECISON THAT **MAKES A MARK**



#### **HANDCRAFTED**

Our highly skilled and experienced surgical mechanics refine and perfect our instruments in Fine Science Tools' very own workshop, a step that only we take in this manner. This is crucial for the precision and consistent quality of our products. For example, our Friedmann-Pearson rongeurs with a tip diameter of 0.5 or 0.7 mm are unique in the world.

Our spring scissors have also established themselves in the research industry and are unsurpassed in their fineness. The result: a level of quality to be proud of – and to which we have now committed ourselves to. Our "Handcrafted" symbol stands for 100% guaranteed craftsmanship and instruments that you can rely on today, tomorrow, and always.






#### **QUALITY CONTROL**

The exceptional quality of our surgical and microsurgical instruments is a direct result of our relentless commitment to detail and craftsmanship. Nearly every instrument we offer is meticulously crafted by skilled European artisans and engineered to exacting specifications.

We use only premium Surgical Stainless Steel and subject each instrument to rigorous testing to ensure outstanding precision, durability, and ergonomic performance. We proudly uphold Supreme German Quality Standards. At our office in Germany, experienced technicians conduct thorough inspections of our instruments.

Every detail is carefully evaluated – from tip sharpness and cutting edges to springs, joints, and overall construction.



#### **CUSTOM DESIGN**

Whether you are seeking a personalized, one-of-a-kind instrument tailored to your specific requirements or interested in building a specialized application kit, we have the expertise and resources to deliver. Through our extensive network of expert partners and collaborations with leading institutions and companies, we offer custom-

manufactured instruments, carefully curated instrument kits, and co-branded sets designed for collaborative projects. Additionally, we provide special modifications to our existing products, ensuring a perfect fit for your unique applications.

Interested? Get in touch with one of our offices - we're happy to assist you!





#### **SCISSORS TEST**

Does paying more for a pair of scissors really buy more value and a longer cutting life?

The picture shows our cutting test for ring scissors based on the standards and specifications of surgical instruments in our office in Heidelberg. Each pair of scissors must perform a smooth cut without tearing or fraying the material, specifically at the fine tips. In addition, a specialized long-term test is done to identify the longevity of other featured alloys and cutting edges, Tungsten Carbide

and ToughCut®. In this test the scissors will repeatedly cut through a continuous feed of cotton cloth and a counter keeps track of the number of cuts made to evaluate when the scissors begin to show a decrease in cutting quality.

The results showcase the addition of added features dramatically increases the cutting life of the scissors. Most notably, the cutting life is more than doubled with the addition of ToughCut® blades and over six times longer with our CeramaCut® feature.

### **OUR RANGE OF SERVICES**



#### **50 YEARS OF EXPERIENCE**



#### **100% SATISFACTION GUARANTEE**



#### PRODUCT AVAILABILITY



#### **FAST DELIVERY TIMES**



#### **TECHNICAL PRODUCT** CONSULTATION





#### **REPAIR AND** MAINTENANCE SERVICE



#### **PERSONALIZED CUSTOMER SUPPORT**



#### **SUPREME GERMAN QUALITY**

STANDARDS



#### TRADE SHOW PRESENCE

IF YOU ARE NOT COMPLETELY SATISFIED WITH YOUR PURCHASE, YOU MAY RETURN IT FOR A FULL REFUND.

> Hans Gawenda, Founder Fine Science Tools Inc.

#### **SHARPENING & REPAIR**

In our very own workshop in Heidelberg, Germany we not only manufacture and customize specific instruments according to customer requirements, but also offer a repair service.

- General repairs
- Sharpening

- Adjustments & corrections
- Cleaning & reconditioning
- Replacements
- Special requests & custom corrections



Our specialized sharpening service offers the highest quality and reliability to restore the original sharpness and functionality of the

Each blade undergoes a thorough inspection for wear, cracks, and other irregularities before sharpening. It is then precisely ground to the optimal angle to ensure superior cutting performance.

With our extensive expertise, we ensure your instruments are repaired with the highest precision, maximizing their longevity and

Each instrument is carefully checked for damage, wear and functionality, followed by a quality check after repair. Our team is also available for repairs of your instruments which have not been purchased from FST®.

For do-it-yourself sharpening, we offer a kit (page 171), which includes a sharpening stone, O-rings, a bottle of instrument oil (50 ml) and a pair of flat nose pliers.



### **ALLOYS & MATERIALS**

There are various special alloys and materials with advanced features, that are used to produce surgical and microsurgical instruments with each having certain properties that would be advantageous in particular

environments or applications. We will introduce these alloys and materials, describe their benefits and highlight unique tools that may be ergonomic and overall helpful for your studies and/or protocols.



#### ToughCut® Scissors

- Identified by black finger rings
- Last approx. 2.3x longer than Surgical Stainless Steel
- One cutting edge has micro serrations to minimize tissue slippage
- The other cutting edge has increased sharpness
- Available as ring and spring scissors



#### **Tungsten Carbide Instruments**

- Identified by short gold-plated handles
- Last approx. 1.9x longer than Surgical Stainless Steel
- Reinforced with tungsten carbide for enhanced durability
- One of the hardest alloys available
- Available instruments: scissors, forceps, needles, needle holders and burrs





#### CeramaCut® Scissors

- Identified by long gold-plated finger rings and a black coating
- Combine Tungsten Carbide and ToughCut® cutting edges with a ceramic coating
- Last approx. 6.6x longer than Surgical Stainless Steel
- Reduced glare
- Reduced surface tension for cleaner cuts
- Increased corrosion resistance



#### Tungsten Carbide & ToughCut® Scissors

- Identified by long gold-plated finger rings
- Last approx. 4.1x longer than Surgical Stainless Steel
- Cutting edges reinforced with tungsten carbide for enhanced durability
- One cutting edge has micro serrations to minimize tissue slippage
- The other cutting edge has increased sharpness



#### Titanium Instruments

- Approx. 40% lighter than Surgical Stainless Steel
- Softer alloy than Surgical Stainless Steel
- Corrosion resistance (cannot rust)
- 100% non-magnetic
- Available instruments: spring scissors, forceps, probes, hooks, spatula and vascular clamps



#### Ceramic coated Instruments

- Identified by full black colour
- Last approx. 2.3x longer than Surgical Stainless Steel
- Increased corrosion resistance
- Reduced glare
- High slide capacity of ceramic reduces surface tension for cleaner cuts
- Available instruments:
   Dumont forceps, scissors



Learn more about our Alloys & Materials





### MAINTENANCE & CARE

Proper care of your surgical instruments will improve their longevity and function. Below are different cleaning, sterilization, and maintenance processes to guide you through your techniques.

- Clean and thoroughly dry your instruments as soon as possible after each use, paying special attention to joints, hinges, and anywhere there is metal to metal contact.
- Ensure that all organic material residues are completely removed.
- Handle your instruments delicately and with care.
- Use your instruments for their intended purpose.
- Store your instruments with tip protectors, and always place them so that they do not touch each other.

#### **CLEANING**

Fine surgical instruments should be cleaned immediately after use. Blood or saline solutions left on instruments for any length of time can quickly start corrosive action. Time, temperature, and agitation play important roles in the cleaning process. Time affects the efficiency of cleaning solutions as most are effective after a certain time period. Higher temperature cleaning solutions result in better cleaning. During manual or ultrasonic cleaning it is helpful to loosen any surgical residue. Before any of the procedures below, rinse your instruments under warm running water to remove all blood and tissue. Instruments can be cleaned manually or ultrasonically, but always:

- Use a neutral pH (pH7) cleaning solution.
- Clean instruments in the open position to ensure a thorough and effective cleaning process.
- Use soft textile or paper cloths, soft plastic brushes made of nylon, or water spray bottle.
- Do NOT use steel wool, wire brushes, or scouring pads as these are abrasive to the instruments.



- Make sure all instruments are visibly clean and free from stains, wiping them carefully to remove all blood and tissue matter.
- · Rinse instruments thoroughly with distilled water, paying close attention to joints and hinges as these areas can accumulate unwanted material.
- Air dry your instruments completely and use a can of compressed air to blow out any water at the joints and hinges. Then store them only in clean and dry environments.



#### **ULTRASONIC CLEANING**

Ultrasonic cleaning is the most effective cleaning method. The vibrating sound waves create micron-sized bubbles in the solution that grow with the alternating pressure. When the bubbles reach a particular size, they collapse creating a force that dislodges dirt and particles, even in the smallest of crevices.



NOT recommended for delicate and micro instruments.

- Follow the manufacturer's guidelines for specific instructions and cleaning cycle times.
- Instruments should be processed for the full recommended cycle time.



- Place instruments in the open position in the ultrasonic cleaner. - Do NOT allow instruments with sharp blades to touch other instruments.
- Do NOT place dissimilar metals (stainless, copper, chrome plated) in the same cleaning cycle.
- Make sure all instruments are completely submerged.
- After the ultrasonic cleaning cycle is complete, rinse instruments thoroughly to remove all ultrasonic cleaning solution.
- Dry instruments thoroughly with a clean towel.
- Use a spray lubricant in the hinges to minimize the risk of corrosion and maintain functionality.



#### STEAM STERILIZATION/AUTOCLAVING

The autoclave process utilizes moist heat (or steam) under pressure for sterilization. The steam used must be from distilled water and free of all impurities.



- Lubricate all instruments that have any metal to metal interaction such as; scissors, rongeurs, needle holders, hemostats, and retractors. Do NOT use WD-40 or other industrial lubricants.
  - Instruments may be autoclaved individually or in sets.
    - Individual instruments: Disposable wraps or sterilization pouches are ideal. Make sure the pouch size is sufficient for the instrument to be in an open (unlocked) position.
    - Instrument sets: Open and unlock all instruments. Place the heavier instruments on the bottom of the set when two layers are needed.



• Do NOT lock instruments during autoclaving. This will prevent the steam from reaching and sterilizing the metal to metal surfaces, and heat expansion during autoclaving can cause cracks in the hinge areas.



- Do NOT over load the autoclave chamber as this may hinder steam penetration.
  - Place a towel on the bottom of the pan to absorb excess moisture during autoclaving. Make sure the towel has no detergent residue and is pH neutral if immersed in water.
  - At the end of the autoclave cycle (before the drying cycle) unlock the autoclave door and open it slightly. Run the dry cycle for the recommended time period by the autoclave manufacturer.

Please note: If the autoclave door is fully opened before the drying cycle, cold air will rush into the chamber, which can cause condensation on the instruments. This will result in water stains on the instruments.



#### **HOT BEAD STERILIZATION**

Recommended for speed and convenience when instrument tips or blade sterilization is required. The Hot Bead Sterilizer does not use dangerous gases, flames, or chemicals. Highlights:

- Bench top set up.
- Only 20 seconds of contact required to ensure your instruments are free of pathogens and
- Beads reach temperature of 240–270 °C (464-518 °F) within 15–20 minutes heating time.

Please see page 174-175 for more information on the FST® Hot Bead Sterilizers.



#### **ETHYLENE OXIDE (ETO) STERILIZATION**

This method is recommended for instruments that will not tolerate high temperatures or moisture (e.g., catheters, cannulae, etc.). This is a process consisting of 4 primary variables; gas concentration, humidity, temperature and time.

Wraps available, see page 176



#### **STORAGE**

Surgical instruments should be properly stored in a clean and dry environment. It is highly recommended to use tip protectors and instrument cases to avoid possible damage to your instruments.

Instrument tip protectors are available, see page 182 Instrument cases are available, see pages 178–181



#### **OXIDATION (RUST)**

All FST® Surgical Stainless Steel instruments are manufactured from medical grade Stainless Steel under strict European standards. It is very important to properly take care of your instruments to extend their lifespan.

Although Surgical Stainless Steel is corrosive resistant, it can still rust or stain if handled improperly. Under unfavorable conditions instruments may oxidize because Surgical Stainless Steel that is exposed to physical, thermal, or chemical manipulations may lose some of its stainless properties.

#### SCISSORS, **BONE INSTRUMENTS, KNIVES & SCALPELS**

Spring Scissors

Fine Scissors	35 – 48	
Surgical Scissors	49 – 55	
Bone Instruments	56 – 69	
Knives & Scalpel	s 70 – 81	

23 - 34



#### **Vannas**

These Vannas-style spring scissors have the smallest blades ever produced. They are ideal where space is restricted and where absolute precision is demanded.



The short taper is ideal for small, precise cuts.

Vannas

MOST Popular





8 cm Cutting edge: 2.5 mm Tip Ø: 0.05 mm

No. 15000-08

Curved

Sharp/Sharp

8 cm

No. 15001-08



8 cm Cutting edge: 2 mm Tip Ø: 0.05 mm

Straight

Sharp/Sharp

No. 15000-03



8 cm Cutting edge: 2 mm Tip Ø: 0.05 mm



Cutting edge: 2.5 mm Tip Ø: 0.05 mm

No. 15000-04



Angled to side Sharp/Sharp

> 8 cm Cutting edge: 2.5 mm Tip Ø: 0.05 mm

> > No. 15002-08

The long taper is ideal for cutting vessel walls or when a longer cutting edge is needed.

Vannas







Straight Sharp/Sharp

8 cm Cutting edge: 3 mm Tip Ø: 0.05 mm

No. 15000-00



Curved Sharp/Sharp

8 cm Cutting edge: 3 mm Tip Ø: 0.05 mm

No. 15000-10

Vannas

Ultra long cutting edge without

sacrificing tip width.

### **SPRING SCISSORS** Vannas-Tübingen MOST Popular Straight Angled up Sharp/Sharp Sharp/Sharp 8.5 cm 8.5 cm Cutting edge: 5 mm Cutting edge: 5 mm Tip Ø: 0.05 mm Tip Ø: 0.05 mm No. 15003-08 No. 15005-08

Angled to side

Sharp/Sharp

8.5 cm

Tip Ø: 0.05 mm

No. 15008-08



**Vannas Highly Polished** 



No. 15403-08


No. 15404-08

**Vannas Round Handled** 



**Vannas with Microserrations** 





**Vannas** 

Curved

Sharp/Sharp

8.5 cm

Cutting edge: 5 mm

Tip Ø: 0.05 mm

No. 15004-08

8.5 cm

No. 15070-08

No. 15071-08

No. 15044-08

No. 15043-08

No. 15009-08

#### **Vannas**

Finished under high magnification, these tips are especially useful when incising vessel walls for catheter insertion.





Moria 9600



Straight Sharp/Sharp

8 cm Cutting edge: 5 mm Tip Ø: 0.15 mm

No. 15396-00



Moria 9601

Curved Sharp/Sharp

8 cm Cutting edge: 5 mm Tip Ø: 0.15 mm

No. 15396-01

No. 15615-08

Straight

Sharp/Sharp

9.5 cm

Cutting edge: 5 mm

Tip Ø: 0.3 mm

No. 15610-08

Curved

Sharp/Sharp


9.5 cm

Cutting edge: 5 mm

Tip Ø: 0.3 mm



· TITANIUM ·



Cohan-Vannas



Straight Sharp/Sharp

6 cm Cutting edge: 6 mm Tip Ø: 0.1 mm

No. 15000-02

Curved Sharp/Sharp

6 cm Cutting edge: 6 mm Tip Ø: 0.1 mm

No. 15000-01



Straight **Blunt/Blunt** 

Angled up

Sharp/Sharp

9.5 cm

Cutting edge: 5 mm

Tip Ø: 0.3 mm

No. 15617-08

6 cm Cutting edge: 6 mm Tip Ø: 0.3 mm

No. 15000-12

Blunt/Blunt 6 cm Cutting edge: 6 mm Tip Ø: 0.3 mm

No. 15000-11



Micro





Sharp/Sharp

10 cm Cutting edge: 7 mm Tip Ø: 0.1 mm

No. 15024-10

Straight

**Blunt/Blunt** 

10 cm

Cutting edge: 7 mm

Tip Ø: 0.3 mm

No. 15025-10



Curved Sharp/Sharp

10 cm Cutting edge: 7 mm Tip Ø: 0.1 mm

No. 15023-10



10 cm Cutting edge: 7 mm Tip Ø: 0.1 mm

No. 15750-11

**Ceramic Coated** 

CERAMIC COATED



Curved Sharp/Sharp

10 cm Cutting edge: 7 mm Tip Ø: 0.1 mm

No. 15751-11



**Ceramic Coated Titanium** 







Sharp/Sharp

11 cm Cutting edge: 5 mm Tip Ø: 0.2 mm

No. 15650-11



11 cm Cutting edge: 5 mm Tip Ø: 0.2 mm

No. 15651-11



Dowell

## a ball tip (Ø 1.3 mm) to prevent piercing of underlying tissue. Angled to side Sharp/Ball tip 10 cm Cutting edge: 8 mm Tip Ø: 0.1 mm Tip Ø: 1.3 mm

Castroviejo with Ball Tip

Same as No. 15006-09 but with











**Vannas** 







Sharp/Sharp 9.5 cm Cutting edge: 2 mm Tip Ø: 0.1 mm No. 15010-09



**Angled Up** 

These spring scissors aid in cutting

tissues from various angles.

and the state of t

No. 15010-10



Tip Ø: 0.1 mm

Straight 11.5 cm



Tip Ø: 0.1 mm No. 15040-11







Angled up Sharp/Sharp

11.5 cm Cutting edge: 10 mm Tip Ø: 0.1 mm

No. 15042-11

No. 15006-09

No. 15033-09

**Vannas Long Handle** 

MOST Popular Straight Curved Sharp/Sharp Sharp/Sharp 10 cm 10 cm

Cutting edge: 4 mm Cutting edge: 4 mm Tip Ø: 0.1 mm Tip Ø: 0.1 mm No. 15019-10 No. 15018-10







Westcott



Pascheff-Wolff


Very delicate blades.

Moria




Moria MC19/B






Noyes











**Noyes Tungsten Carbide** 

Tungsten carbide cutting edges provide long lasting sharpness.

TUNGSTEN



12 cm Cutting edge: 14 mm Tip Ø: 0.2 mm

No. 15514-12

No. 15011-12

Tip Ø: 0.2 mm

No. 15012-12

Micro 2000 ToughCut®





Straight Sharp/Sharp

12.5 cm Cutting edge: 6 mm Tip Ø: 0.3 mm

No. 15124-12



Curved Sharp/Sharp

12.5 cm Cutting edge: 6 mm Tip Ø: 0.3 mm

No. 15123-12 No. 15400-12

Straight

Sharp/Sharp

12.5 cm

Cutting edge: 6 mm

Tip Ø: 0.05 mm

#### **Vannas Round Handled**

Round handles allow easy fingertip adjustments.











12.5 cm Cutting edge: 6 mm Tip Ø: 0.05 mm No. 15401-12



14 cm Cutting edge: 6 mm Tip Ø: 0.1 mm

No. 15020-15

14 cm

Cutting edge: 6 mm

Tip Ø: 0.1 mm

No. 15021-15

For microscopic dissections of soft and mucous tissue. The minute fineness of the blades allows clear visualization of the tissue being dissected.

Clipper



#### Uses:

- Removing sheath of ganglia to access neurons in aplysia.
- Cutting single nerve axons in aplysia, mollusk and snail.
- Stripping fascia off arteries and muscles.
- Ultra fine brain tissue dissection.
- · Cutting of minute fiber optic strands.



8 cm Cutting edge: 3 mm Tip Ø: 0.05 mm

Blade shape: Triangular

No. 15200-00

#### Trident II

Special design with short handle and inner spring.



Straight Sharp/Sharp

8 cm Cutting edge: 8 mm Tip Ø: 0.2 mm



8 cm Cutting edge: 8 mm Tip Ø: 0.2 mm

No. 15110-08

No. 15111-08



PRICES SUBJECT TO CHANGE PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER

Micro with Pen-Type Handle

The finest blades on spring scissors. Designed for ophthalmic With microserrations to minimize slippage.

### surgery (vitrectomy). The spring handles are made of ribbed plastic to ensure a solid grip. Autoclavable.

Micro with Pen-Type Handle





Bonn-Gawenda

The innovative Bonn-Gawenda Scissors from our Exclusive Line combine a precise cutting edge,

Bonn

Our finest iris scissors, best suited for microdissection under high magnification.







9 cm Cutting edge: 10 mm

No. 14084-08



9 cm Cutting edge: 10 mm



9 cm Cutting edge: 10 mm

No. 14085-08 No. 14083-08

#### F·S·T

Straight

Sharp/Sharp

14.5 cm

Cutting edge: 2.5 mm

Tip Ø: 0.3 mm

No. 15030-14

Angled to side

Sharp/Sharp

14.5 cm

Cutting edge: 2.5 mm

Tip Ø: 0.3 mm

No. 15031-14

**Bonn Highly Polished** 

# Highly polished. Straight Curved Sharp/Sharp Sharp/Sharp 9 cm 9 cm





**Bonn-Strabismus ToughCut®** 








Iris Sandblasted

and Vanadium.





No. 14061-11

Iris

Delicate pattern.

No. 14085-09



**Moria Iris** 










**Iris Large Loops** 





No. 14041-10



Sharp/Sharp

12 cm

Cutting edge: 20 mm

No. 14069-12

Sharp/Sharp


12 cm

Cutting edge: 20 mm

No. 14068-12



Wagner



12 cm Cutting edge: 20 mm

No. 14070-12 No. 14071-12

Cutting edge: 20 mm

F·S·T

No. 14160-10

No. 14370-22

No. 14161-10

No. 14040-10

### Iris ToughCut® TOUGH MOST Popular Straight Curved Sharp/Sharp Sharp/Sharp 9 cm 9 cm Cutting edge: 16 mm Cutting edge: 16 mm No. 14059-09 No. 14058-09 11.5 cm 11.5 cm Cutting edge: 16 mm Cutting edge: 16 mm No. 14058-11 No. 14059-11





**Iris Tungsten Carbide** 

**POPULAR** 

TUNGSTEN

CARBIDE



Iris ToughCut® and Tungsten

Carbide

MOST Popular









Iris

Straight Curved



Lexer-Baby








Metzenbaum Mini

TUNGSTEN



No. 14574-09



**Metzenbaum Baby Tungsten Carbide** 

TUNGSTEN

CARBIDE



11 cm Cutting edge: 20 mm Cutting edge: 20 mm No. 14574-11

No. 14575-11

Sharp/Blunt

10 cm

Cutting edge: 17 mm

No. 14028-10

Sharp/Blunt

10 cm

Cutting edge: 17 mm

No. 14029-10

No. 14575-09

Shea

Transplant scissors with light and delicate pattern.





12 cm Cutting edge: 16 mm

No. 14105-12

Cutting edge: 20 mm

No. 14074-11

Cutting edge: 10 mm

No. 14077-09

Cutting edge: 10 mm

No. 14077-10

Cutting edge: 20 mm

No. 14075-11

Cutting edge: 13 mm

No. 14067-11

Cutting edge: 13 mm

No. 14066-11

Cutting edge: 12 mm

No. 14173-12

Similar to 14081-09, but with stronger blades for cutting fine bone.

Used for cutting stitches or suture.

Iris scissors with one ball tip to prevent piercing of underlying tissue.



**FINE SCISSORS** 















Angled to side Sharp/Sharp

10 cm Cutting edge: 6 mm

No. 14081-09



Angled to side Sharp/Sharp

10 cm Cutting edge: 8 mm

No. 14082-09

Straight Blunt/Blunt

9 cm Cutting edge: 16 mm

No. 14076-09



Straight Sharp/Ball tip

9 cm Cutting edge: 10 mm 1.5 mm Ball Diameter

No. 14086-09



Straight Sharp/Ball tip

11.5 cm Cutting edge: 16 mm 1.5 mm Ball Diameter

No. 14080-11



Straight **Blunt/Blunt** 

11.5 cm Cutting edge: 16 mm

No. 14072-10



Angled up Sharp/Sharp

10 cm Cutting edge: 12 mm

No. 14037-10

# Stevens with Adjustable Loops LARGE LOOP Straight Blunt/Blunt 10 cm Cutting edge: 8 mm No. 14108-09

#### Iris with Adjustable Loops The finger loops can be enlarged to fit almost any size.

LARGE LOOP

#### Kaye with Adjustable Loops

These scissors have semi-sharp outer edges to better facilitate dissection.









10.5 cm Cutting edge: 14 mm

No. 14171-10

#### **Standard Surgical**

Our standard surgical scissors are recognized for these important features:

- Excellent feel and handling capability.
- Overall execution that marks them as fine examples of German craftsmanship.





	14008-14				
	Length (cm)	Cutting edge (mm)	Straight Sharp / Sharp	Straight Sharp / Blunt	Straight Blunt / Blunt
	12	26	No. 14002-12	No. 14001-12	No. 14000-1
	13	28	No. 14002-13	No. 14001-13	No. 14000-1
	14.5	30	No. 14002-14	No. 14001-14	No. 14000-1
	16.5	36	No. 14002-16	No. 14001-16	No. 14000-1
	18.5	40		No. 14001-18	No. 14000-1
0	20	44			No. 14000-2
2	Length (cm)	Cutting edge		Straight Sharp / Blunt	
F-S-T 14	14.5	30		No. 14008-14	
T 14					

Length (cm)	Cutting edge	Curved Sharp / Sharp	Curved Sharp / Blunt	Curved Blunt / Blunt
12	26	No. 14005-12	No. 14004-12	No. 14003-12
13	28	No. 14005-13	No. 14004-13	No. 14003-13
14.5	30	No. 14005-14	No. 14004-14	No. 14003-14
16.5	36	No. 14005-16	No. 14004-16	No. 14003-16
18.5	40		No. 14004-18	No. 14003-18
20	44			No. 14003-20

Sharp/Sharp 9.5 cm Cutting edge: 8 mm

No. 14106-09

Lister

Excellent skin cutting scissors. Longer, blunt blade prevents damaging tissue underneath skin without losing the cutting line.







Blunt/Blunt

14 cm Cutting edge: 26 mm

No. 14131-14

**Metzenbaum Fino** 

Metzenbaum ToughCut®

TUNGSTEN CARBIDE 70% of 80% of actual size actual size Straight Curved **Blunt/Blunt** Blunt/Blunt 13 cm 13 cm Cutting edge: 20 mm Cutting edge: 20 mm Straight Straight Curved Curved No. 14018-13 No. 14019-13 **Ultra fine Blunt/** Blunt/Blunt **Ultra fine Blunt/** Blunt/Blunt Blunt **Blunt** 14.5 cm 14.5 cm 14.5 cm 14.5 cm 14.5 cm 14.5 cm Cutting edge: 22 mm No. 14518-14 No. 14018-14 No. 14019-14 No. 14508-15 No. 14509-15 No. 14519-14 18 cm 18 cm 18 cm 18 cm Cutting edge: 24 mm Cutting edge: 24 mm Cutting edge: 24 mm No. 14018-18 No. 14019-18 No. 14518-18 No. 14519-18 PRICES SUBJECT TO CHANGE PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER

**Metzenbaum Fino Tungsten Carbide** 



Metzenbaum




#### Friedmann-Pearson

Experience the unparalleled precision of our flagship rongeur which features a specialized 0.5 mm cup diameter, tailored for the finest research techniques. These rongeurs guarantee effortless cutting while ensuring captured tissue remains secure and intact.









#### Friedmann-Pearson

Best suited for delicate work, this is one of our finest tip rongeurs in our product line. The fine-cupped jaws have been designed with a longer and slimmer approach for better visibility.







14 cm

Cup size: 0.7 mm

No. 16120-14

Curved 14 cm

Cup size: 0.7 mm

No. 16121-14

14 cm Cup size: 1 mm No. 16020-14

Straight




14 cm Cup size: 1 mm

No. 16021-14

Pearson

Provides excellent visibility at dissection site due to its extra long and narrow approach to the fine cutting jaws. Especially suitable for larger hands.





Luer





Curved


12.5 cm

Cup size: 2.5 mm

No. 16000-14

Extra fine jaws.

Lempert



**Shallow bend** 

17 cm

Cup size: 1 mm

No. 16015-17

Sharp bend

17 cm

Cup size: 1.2 mm

No. 16016-17

Fine

llic

Round cutting blades for tissue, cartilage and soft bone.



Cutting edge on both outside and inside bevel. Same structure and functionality as Item No. Cuts suture, cartilage, bone and steel wire to 11215-01 (left), but with an epoxy coating up to 0.06 mm, and copper wire to 0.15 mm. the cutting edge for easier care and maintenance.





### F·S·T by DUMONT 9

**Dumont #15ARW Cutter** 

### NOT AUTOCLAVABLE



Tungsten Carbide #15AGHM Cutter

High precision cutting tweezers are suitable for applica-

tions which require cutting small material or performing

small cuts. The extreme hard Tungsten Carbide blade







11.5 cm Cutting edge: 7 mm Jaw Opening: 2.5 mm

Epoxy coated carbon steel

No. 11215-02



**AngledCutting** 11.5 cm

Cutting edge: 10 mm Jaw Opening: 4 mm Anti-magnetic anti-acid Surgical Stainless Steel (Inox) (sa), Tungsten Carbide, sandblasted

No. 11515-02





14 cm Cutting edge: 8 mm

No. 16109-14





Curved

14 cm Cutting edge: 9 mm

No. 16107-14



**AngledCutting** 

11.5 cm

Cutting edge: 10 mm

Jaw Opening: 2.5 mm

Carbon steel

No. 11215-01

**Bone Nipper** 

Cohen

Strong jaws for cutting tough cartilage.







10 cm Cutting edge: 7 mm

No. 16101-10









11 cm Cutting edge: 9 mm No. 16140-11



Straight

Blunt/Blunt


13 cm

Cutting edge: 12 mm

No. 16144-13

**Spring Scissors Style with** Strong Handle

Slim, strong blades for cutting fine bones and tissues.





Curved Blunt/Blunt

12 cm Cutting edge: 12 mm

No. 16143-12



11 cm Cutting edge: 5 mm

No. 16102-11



Littauer

Straight 11 cm Cutting edge: 15 mm No. 16152-11 Straight Straight 12.5 cm 15 cm



Beebee

Very strong jaws with sharp points.



No. 16030-16

No. 16030-12

**Scissors Tungsten Carbide** 

Strong heavy blades for cutting bone.



Boehler Bone splitting cutter with very fine tips and double

action joint to relieve strain on your hand.



**Skull Flap Cutter** 

One tip has a flat and narrow platform which fits be-


tween dura and skull. Other tip has sharp cutting edge. The double action joint relieves strain on your hand.

Cleveland

Gross anatomy bone cutter.

Liston

**Bone Pliers** To peel back skull flaps. **Bone Shears** 



17 cm

Cutting edge: 15 mm


No. 16105-16



20 cm

Cutting edge: 30 mm

No. 16104-19





Curved

24 cm Cutting edge: 65 mm

No. 16150-24

### **Bone Forceps**

Strong cross-serrated tips for bone crushing and removal.



# Angled

11.5 cm Tip Ø: 1 mm

No. 16060-11

### **Mini Self-Tapping Bone Screws**

Used extensively to anchor dental cement to bone, these self-tapping screws have superior holding strength over ordinary machine screws. Surgical Stainless Steel with fillister flat head.



Shaft diameter: 0.86 mm Pkg. of 100 screws

No. 19010-10

Pkg. of 10 screws

No. 19010-11

For pre-drilling use our 0.7 mm burrs, see page 183.

Shaft diameter: 1.19 mm Pkg. of 100 screws

4.8 mm

No. 19010-00

Pkg. of 10 screws

No. 19010-01

For pre-drilling use our 0.9 mm burrs, see page 183.


### **Bone Wax**


A beeswax-based compound for bone hemostasis. Soft, kneadable and ready to use. Sterile.



Pkg. of 12 foil packs (2.5 g each)


No. 19009-00





### **Gross Anatomy Knives**

Handles are made of fiberglass reinforced yellow nylon with a slip resistant surface that allows little chance of bacterial accumulation. The blades are made of German Molybdenum Stainless Steel.



### **Dissection knives**

26 cm Cutting edge: 13 cm

No. 10157-13

28 cm Cutting edge: 15 cm

No. 10157-15

31 cm

Cutting edge: 30 cm No. 10152-30

44 cm

**Brain sectioning knife** 

40% of

actual size

### **Obsidian Stone Age Scalpels**

For applications where extremely fine cuts are required or for studies where trace metals from ordinary scalpel blades cannot be tolerated. These ultra fine cutting edges have been used in human surgical procedures where objectionable scar tissue results from the use of standard scalpels. The blade is made from obsidian — a type of volcanic glass — which allows a much finer blade than is possible with the conventional steel blade. In some cases this makes the relatively inexpensive obsidian scalpel a suitable substitute for an expensive diamond knife. Patterned after "stone age" knives which first made their appearance about 18,000 years ago, these modern versions utilize hard wood and plastic dowels for handles. Plastic dowels are made of polyoxymethylene.



Obsidian is a very fragile material. Great care should be taken to not exert any lateral pressure on the blade during cutting. Each blade is hand-fashioned, so sizes, shapes and points will vary. While all the blades will have at least one sharp edge, some will have two.



Cutting edge: 18 cm

No. 10157-18

Straight 13 cm

Cutting edge: 6 mm No. 10110-01

Cutting edge: 12 mm No. 10110-03

13.5 cm

12.5 cm Cutting edge: 6 mm

No. 10110-05

12.5 cm Cutting edge: 6 mm

No. 10110-07

### **Ceramic Blades and Handle**

- Ceramic edge stays sharp 100 times longer than steel edge
- Rust and corrosion free
- Chemical and acid resistant
- Withstands up to 2400 °C (4352 °F)
- No electrostatic discharge concerns
- No trace metal contamination
- Fully nonmagnetic



Advantages:

This knife has a gem-quality, neutral diamond blade which is honed to the finest cutting edge possible. The blade is made to an ophthalmic surgical standard and is inspected under 400 x magnification. The retractable handle aids in care and storage.

**Diamond Knives** 



This blade is extremely fragile and vulnerable to lateral pressure.





### Angled 45°

33 x 6 x 0.5 mm Cutting edge: 4.5 mm Pkg. of 5 blades

No. 10025-44

### Angled 45°

33 x 6 x 0.5 mm Cutting edge: 2.5 mm Pkg. of 5 blades

No. 10025-43

### Angled 23°

33 x 6 x 0.5 mm Cutting edge: 14 mm Pkg. of 5 blades

No. 10025-23

### Angled 45°

12.5 cm

No. 10025-50




33 x 6 x 0.5 mm

No. 10025-45









## Lancet Cutting angle: 30°

12 cm Cutting edge: 1 mm Tip thickness: 0.15 mm


No. 10100-30

### Angled Cutting angle: 45°

12 cm Cutting edge: 1 mm Tip thickness: 0.15 mm

No. 10100-45

### **Blade Breakers and Holders**



No. 10052-11



Concave / Convex surface Handle diameter: 5.5 mm

13 cm Jaw width: 4 mm Concave / Convex surface

Handle diameter: 8 mm No. 10053-13

### **Breakable Blades**

These special, high carbon content blades are the only blades that will break cleanly.






Blade thickness: 0.1 mm Pkg. of 10 blades

No. 10050-00


### **Create Your Own Blade**

With these blade holders/breakers and our special carbon steel blades, No. 10050-00, you can produce microblades in sizes that meet your own individual needs.



Aluminum Alloved Handle

Cutting edge: 8 mm Pkg. of 5 blades





No. 10053-09

KNIVES & SCALPELS

**SCALPELS** 

### **Sterile Scalpels**

Cold forged from carbon steel and sharpened to eliminate tissue distortion and blade drag. 12.5 cm plastic handle with slip-proof grip. Sterile.









**Crescent blade** 

14 cm Cutting edge: 11 mm Tip thickness: 0.2 mm

No. 10317-14



15° Cutting Angle

13 cm Cutting edge: 5 mm

Tip thickness: 0.15 mm

No. 10315-12

4.0 mm

Needle blade

14.5 cm

Cutting edge: 4 mm

Tip thickness: 0.2 mm

No. 10318-14

Handle for Inter-

changeable Blades 9.5 cm

No. 10035-00

12.5 cm

No. 10035-12



**Interchangeable Blades and Handles** 



Angled Carbon steel

Cutting edge: 10 mm Pkg. of 5 blades

No. 10035-15



Curved

Carbon steel

Cutting edge: 10 mm Pkg. of 5 blades

No. 10035-20

Small / Flat Carbon steel

Cutting edge: 1.5 mm Pkg. of 5 blades

Cutting edge: 2.75 mm Pkg. of 5 blades

No. 10035-10

Large / flat Carbon steel

No. 10035-05



Cane

Carbon steel

Cutting edge: 12 mm Pkg. of 5 blades

No. 10035-25

Arc Carbon steel

No. 10035-30

### Standard Scalpel Handle #4 with Silicone

For heavy, gross anatomy work. Molded to fit



Especially suitable for large dissections.

The alloy used has a high carbon content which provides a longer lasting edge. Blades come individually wrapped. Non-sterile.

**Blades for #4 Scalpel Handles** 





#20 blades

Box of 100 blades

No. 10020-00



#21 blades

Box of 100 blades

No. 10021-00



#22 blades

Box of 100 blades

No. 10022-00



#23 blades

Box of 100 blades

No. 10023-00



Cutting edge: 12 mm Pkg. of 5 blades

No. 10009-13

16 cm

13 cm

No. 10004-13

KNIVES & SCALPELS

### Standard Scalpel Handle #3

MOST Popular

### Secure Scalpel Handle #3

1 Slide blade onto front pin of the handle.

2 Then slide pin to front slot

of blade. Rotate blade clockwise and secure blade to back pin on handle.

**3** Scalpel is ready to use.

A safety lock holds the blade securely, and allows The alloy used has a high carbon content which simple and safe blade attachment and removal. provides a longer lasting edge. Blades come individually wrapped. Non-sterile.











Blades for #3 and #7 Scalpel Handles

### #10 blades Box of 100 blades

No. 10010-00



Box of 100 blades

No. 10011-00



### #12 blades


Box of 100 blades

No. 10012-00



### #15 blades

Box of 100 blades



**Short Scalpel** Handle #7

**Standard Scalpel** 

Handle #7

### **Scalpel Blade Remover**


For safe removal of scalpel blades, simply insert scalpel with blade into the Scalpel Blade Remover and press the sides together. This will dislodge the blade from handle and hold the blade securely for release.



### **Scalpel Blade Remover and Disposal Box**

This plastic device will remove used scalpel blades safely and retain them for safe disposal. Can be used for blades that fit the #3, #4, #7 and Gross Anatomy handles.







F.S.T 10003-13

F·5·T 10003-12



12 cm No. 10003-12

F·S·T

No. 10003-13

12 cm



No. 10015-00



12 cm

No. 10007-12

16.5 cm

No. 10007-16

KNIVES & SCALPELS

### **Disposable Safety Scalpels**

Plastic handle #3 with Surgical Stainless Steel blades. The protective cartridge completely covers the blade for safe, fast and easy disposal. Single-use and non-sterile.





Scalpel shown with safety cartridge closed and locked into place.

### #60 blades 7.5 cm 13 cm

Box of 100 blades No. 10060-00

#10 scalpels 14 cm Pkg. of 10 blades

No. 10000-10


#11 scalpels 14 cm Pkg. of 10 blades No. 10000-11

14 cm Pkg. of 10 blades No. 10000-15

#15 scalpels

### **Safety Scalpel Blade Cartridges**

The Safety Scalpel Handle #3 is compatible with all Safety Scalpel Blades. The handle is reusable and made of Teflon coated stainless steel which provides a non-stick surface. The Blade Cartridges are color-coded with clear windows for easy identification. The protective cartridge completely covers the blade for safe, fast, and easy disposal. Handle is reusable and autoclavable. Blades are single-use and non-sterile.





12.5 cm

Teflon coated Surgical Stainless Steel (Inox)

No. 10006-12

sold separately.

#10 blades

Pkg. of 10 blades

No. 10010-10



#11 blades

Pkg. of 10 blades

No. 10011-10

#15 blades Pkg. of 10 blades

No. 10015-10

F·S·T

No. 10060-13

F-5-T

**FST® BY DUMONT** 

**FORCEPS, TWEEZERS, HEMOSTATS** 

83 - 92

97 - 98

98, 100

104 - 115

116 - 121

122 - 131

93 - 96, 99 - 103

& VASCULAR INSTRUMENTS

FST® by Dumont

Fine

Moria

S&T®

Standard

Hemostats

Vascular Instruments


These forceps have been honed

### **Dumont #5SF**

The finest tips ever produced by Dumont.

MOST POPULAR

# Finished under magnification and



**Dumont #5 Medical Biology** 





**Dumont #5 Ceramic Coated** 

These forceps are non-reflective.

The scratch and corrosion resistant

coating guarantees durability.

CERAMIC COATED





### 11 cm 0.025 x 0.005 mm Surgical Stainless Steel (Inox)

No. 11252-00



11 cm 0.05 x 0.02 mm Surgical Stainless Steel (Inox)

No. 11254-20



11 cm 0.05 x 0.02 mm Ceramic coated Surgical Stainless Steel (Inox)

No. 11252-50




### **Biology**

11 cm 0.05 x 0.02 mm Dumostar

No. 11295-20





Straight 11 cm 0.1 x 0.06 mm

Steel (Inox) No. 11253-20

Surgical Stainless

Standard Off-Center 11.5 cm 0.13 x 0.08 mm Surgical Stainless

Steel (Inox) No. 11253-21






Curved 11.5 cm 0.17 x 0.1 mm Surgical Stainless

Steel (Inox) No. 11273-22






Curved 11.5 cm 0.17 x 0.1 mm Surgical Stainless Steel (Inox)

No. 11273-20


Medical #7

Medical #5/45



11 cm 0.1 x 0.06 mm Surgical Stainless Steel (Inox) No. 11253-25

Medical #5/90



11 cm 0.1 x 0.06 mm Surgical Stainless Steel (Inox)

No. 11253-29

### Dumont #55

F·S·T

by DUMONT 9

Bevelled, like #5 but with extra light and fine shanks.

### **Dumont #5 Mirror Finish**

The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves or any

other adhesive tissue.






55

11 cm 0.05 x 0.02 mm Surgical Stainless Steel (Inox)


# No. 11255-20

11 cm 0.05 x 0.02 mm Dumostar

No. 11295-51



11 cm 0.1 x 0.06 mm Surgical Stainless Steel (Inox) No. 11251-23



0.05 x 0.02 mm Surgical Stainless Steel (Inox) No. 11252-23

### Dumont #5

These are the most widely used Dumont forceps. Bevelled, very fine shanks for best visibility.



0.10 mm

Standard

11 cm

0.1 x 0.06 mm

Surgical Stainless

Steel (Inox)

No. 11251-20

Dumoxel

No. 11251-30


Dumostar

No. 11295-00















Dumostar

No. 11295-10 11 cm 0.08 x 0.04 mm

0.08 mm

Standard

Carbon steel

No. 11251-10


Titanium No. 11252-40

**Dumont #5/45** 

is squared for a better grip.

**Dumont #5XL** 





No. 11251-21

<u>0.05 mm</u>


**Biology** 

12 cm

No. 11252-21









F·S·T

by DUMONT 9

Dumont #L5



Dumont #5

Teflon® Coated

These antimagnetic forceps

are ideal for use in cryogenic

**Dumont #5/45** 

 $45^{\circ}$  tips for improved visibility and

control.



No. 11251-35

No. 11626-11

No. 11251-33

F·S·T

by DUMONT 9

No. 11231-30

No. 11293-00

No. 11274-20

No. 11297-00

No. 11272-40

No. 11270-20

No. 11272-50

No. 11260-20

No. 11231-20

**Dumont SS** 

specimen is ensured.

### Dumont #2

These forceps are ideal for breaking off delicate bones.

### **Dumont #2AP Epoxy Coated**


Provides good insulation for up to 220 AC voltage work and highly impervious to chemical solutions. Frozen or super cold sections will not adhere to the gripping surfaces.


















F·S·T

by DUMONT 9

**Dumont AA** 



**Dumont AA Epoxy Coated** 

Slip-proof acid and electric shock

resistant epoxy coating covers the

forceps 17.5 mm from the tips.



0.13 mm

Standard 11 cm

0.13 x 0.08 mm

Dumoxel

No. 11241-30

Dumostar

No. 11294-00

0.06 mm

**Biology** 

11 cm

0.06 x 0.04 mm

Dumoxel

No. 11242-40

**Dumont Mini Forceps** 

**Dumont #5AC** 

Balanced #5

These forceps feature an ergonomic design and ultra fine tips for optimal comfort and accuracy. The balanced action of the forceps ensures seamless efficiency during prolonged tasks, providing ease in rotational pivoting with the rounded handle and a counterweight technology providing enhanced control and stability.

**Gold Plated** 


Gold-plating gives corrosion resistance to most chemicals, salts and acids. These forceps are perfect for any application where corrosion resistance and toughness are primary requirements.















0.12 mm Standard 7 cm 0.12 x 0.07 mm


Mini #3

7 cm Surgical Stainless Steel (Inox)

\_\_\_\_0.14 mm

Mini #4

Standard 0.14 x 0.06 mm Surgical Stainless Steel (Inox) No. 11200-12



Mini #5

Surgical Stainless Steel (Inox) No. 11200-14

0.1 x 0.06 mm Gold plated

Surgical Stainless Steel (Inox)

No. 11915-11

Gold plated

Surgical Stainless Steel (Inox) No. 11917-12

Steel (Inox) No. 11911-11

11 cm

Surgical Stainless

No. 11200-10

### **Diamond Coated**

Diamond coating provides increased corrosion resistance and greater protection from wear and tear, which extends the operational life of the forceps, as well as

providing a stronger grip at the tips.

### #27 Mirror Finish

These forceps have long and narrow shanks with thin, highly polished tips to prevent tissue from sticking. Ideal for embryos, nerves, or other adhesive tissues.

### #259 Replaceable Tips


Surgical Stainless Steel body with replaceable plastic tips. Tips are made of polyvinylidene fluoride (PVDF) and are electrostatic discharge (ESD) safe. Body and plastic tips are corrosion and abrasion resistant.

### #5 Mirror Finish

The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves or other adhesive tissue.

### #3

Straight forceps with sturdy coarse tips.

















11 cm 0.1 x 0.06 mm 0.05 x 0.02 mm Surgical Stainless


Surgical Stainless Steel (Inox)

Steel (Inox) No. 11412-11 No. 11413-11

0.05 mm


Straight

11 cm



12 cm 0.2 x 0.12 mm Surgical Stainless Steel (Inox)

No. 11445-12



0.15 x 0.1 mm Diamond coated Surgical Stainless Steel (Inox)

No. 11990-14

0.15 mm



0.15 mm






Moria Iris

Balanced with gentle spring action and tips

that are perfectly matched.









F·S·T 11992-20








### Moria #5

These forceps have the most finely executed points available. They are extremely light, superbly balanced and have a gentle spring action. The shanks are specially designed with an arch which ensures maximum force to the tips. The tips close flat for 6 mm from the points to give a large holding area. The tips have been expertly finished under magnification.









**Straight** 

13 cm

0.1 x 0.06 mm

No. 11399-80


NOT FOR HUMAN USE PRODUCTS ARE NOT INTENDED FOR USE ON OR IN HUMANS - MOST ITEMS SHOWN IN ACTUAL SIZE

Curved

13 cm

0.1 x 0.06 mm

No. 11399-87



# **Moria Paufique** Moria Moria 3225


Magnified View

### **#5 Fine Forceps**

S&T®

The tips are processed with MicroCrystals<sup>®</sup>, a method developed by S&T® to infuse Tungsten Carbide on Surgical Stainless Steel. The formation of this ultra-fine coating guarantees atraumatic handling of tissue while improving the grip dramatically. The recommended use is where forceps are utilized as needle holders or when

#5 SuperGrip Forceps



### Bonn **Castroviejo Suture Tying**

The 1x2 teeth forceps No. 11084-07 are the finest

of their kind with 0.1 mm wide teeth. The smooth


forceps No. 11083-07 have a narrow 5 mm long

platform. Both forceps have a pin inside the

shanks to limit the applied pressure.

Magnified View

The fine Bonn forceps features a tying platform and is also available with additional 1x2 teeth. It's designed for tissue manipulation and suturing. The wide handle ensures safe handling. The fine and precise tips make it suitable for work under a high magnification.




### **Mirror Finish Suture Tying Forceps** with Teeth

A unique tip style fusing a suture tying platform and a 1x2 teeth tip for tissue manipulation.











S&T JF-5





S&T JF-5 TC






No. 11084-07







No. 11090-09

No. 11390-25

No. 11083-07

### **Moria Troutman Suture Tying**

Used to suture after surgery is complete. This requires the manipulation of a calibrated nylon thread and the use of a specific instrument: a tying forceps, straight or curved, without teeth.









Straight Smooth with platform

9.5 cm Tip Ø: 0.3 mm Tying platform: 8 mm

No. 11393-04

### **McPherson**

MOST Popular

0.4 mm

Angled 45°

Smooth with

platform

9 cm

Tip Ø: 0.4 mm

Tying platform: 6 mm

No. 11063-07

The platform tips close flat. This ensures a reliable grip on sutures.


**#3 Suture Tying Forceps** 



# **Tuebingen Suture Tying**

The platform tips close flat. This ensures a reliable grip







13.5 cm


Tip Ø: 0.3 mm

Tying platform: 6 mm

No. 00272-13

# on sutures.







### **Round Handled Suture Tying**

Round handles allow easy fingertip adjustments.





Smooth with platform 10.5 cm

Tying platform: 5.5 mm

No. 18026-10

### **Vial Forceps**

Tips taper for 6 cm allowing easy access to specimens in vials.





15 cm Tip Ø: 0.5 mm

No. 11006-15

F·S·T

PRICES SUBJECT TO CHANGE PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER

Moria 9605

Curved

Smooth with

platform

8.5 cm


Tip Ø: 0.3 mm

Tying platform: 7.5 mm

No. 11393-05

**Graefe Forceps** 


Graefe



Bishop-Harman

Foerster

Fine Graefe



### **Ring Forceps**

For transport of tissue and organs. With limit screw to prevent excessive applied pressure.



atraumatic gripping of

organs and tissues.

### **Live Insect Handling**

Very delicate.

### Micro-Adson

Our most delicate version of the popular Adson forceps.

### **Micro-Adson with Fenestrated Handle**

Our delicate Micro-Adson forceps, with five windows to aid in stability and easy finger adjustments.

### Adson

Wide, strong grip with narrow tips.











Straight

Serrated

12 cm



Straight

12 cm

1.5 x 2 mm

Straight 1x2 Teeth **Cross Serrated** 

F-5-T 11006-12

**Tungsten Carbide** 12 cm 1.5 x 1.3 mm


No. 11506-12















►I | 0.8 mm

F.S.T 11018-12



12 cm 0.5 x 0.6 mm

Straight

Serrated

12 cm

0.8 x 1 mm

No. 11020-12

F-S-T

No. 11022-12

Ring Ring Ring 8.5 cm 9 cm 9 cm ID: 1 mm OD: 1.5 mm ID: 2.2 mm OD: 3 mm ID: 4.8 mm OD: 6 mm x x 0.5 mm x 1.3 mm 1.5 mm

No. 11103-09

Straight



No. 11072-11



Straight

No. 26029-10



1x2 teeth 12 cm 0.8 x 2 mm Tying platform: 6 mm

Straight

Tying platform with

No. 11027-12

Straight

No. 11101-09

Straight

No. 11106-09

### Adson DeBakey

Ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures. Two rows of atraumatic teeth interfit with one opposite row of atraumatic teeth.








### Atraumatic **DeBakey Serrations**

This tooth pattern, as shown, consists of one row of atraumatic teeth meshed with two opposite rows of atraumatic teeth.

These two forceps are ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures.




Adson-Graefe



4.5 x 1.25 mm No. 11030-12



Adson-Brown



No. 11627-12


### **Anatomical Standard Models**

Broad pattern with strong spring action.



Straight Serrated

ILEIII	(mm)	(cm)
No. 11000-12	2.5	12
No. 11000-13	2.8	13
No. 11000-14	3.25	14.5
No. 11000-16	3.3	16
No. 11000-18	3.7	18
No. 11000-20	3.75	20
No. 11000-25	3.9	25
No. 11000-30	4	30



Length (cm)	lip width (mm)	Item
12	2.5	No. 11001-12
13	2.8	No. 11001-13
14.5	3.25	No. 11001-14
16	3.3	No. 11001-16
18	3.7	No. 11001-18
20	3.75	No. 11001-20

### **Anatomical Narrow Models**

Lighter spring action than standard pattern for grasping more delicate tissue.





Straight

11602-14



F-S-T 11002-13



3.2



No. 11002-25

Length (cm)	Tip width	Item
12	2	No. 11003
13	2.15	No. 11003
14.5	2.25	No. 11003
16	2.35	No. 11003
18	2.4	No. 11003
20	2.75	No. 11003

F·S·T

Straight

DeBakey

serrations

12.5 cm

1.75 x 1.5 mm

No. 11617-12

Straight

DeBakey

serrations

15 cm

1.75 x 1.5 mm

No. 11616-15

F-S-T 11000-13

### Semken

The Semken tweezers are characterized by their slim shape and fine tips. Their anatomical (toothing) or surgical (1x2 teeth) tips make them ideal for grasping tissue. The serrated grip surface ensures a secure hold during use.





Straight Serrated 13 cm

Tip width: 1.3 mm No. 11008-13

15 cm Tip width: 1.5 mm

No. 11008-15

15 cm Tip width: 1.5 mm Titanium

No. 11608-15

### **DeBakey Round Handled**

Ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures. Two rows of atraumatic teeth interfit with one opposite row of atraumatic teeth. The round handles allow easy fingertip adjustments.



Straight **DeBakey serrations** 

15 cm Tip width: 1.2 mm

No. 11076-15

### **Multipurpose Forceps**

Balanced instruments are designed to put the center of gravity between the thumb and index finger. This is accomplished by adding a slight counter weight to the proximal end of the instruments. The design enhances tactile feel, aiding in better control and precision.


The tips of these multipurpose forceps are configured to enable their use as dilators at the front, as well as needle holders and tying forceps further in the back. The platform at the tips close parallel ensuring a reliable grip.



Straight For Dilating and Suture tying

14.5 cm Tip Ø: 0.1 mm

No. 18601-15



Straight

For Dilating and Suture tying

15 cm

Tip Ø: 0.3 mm

No. 18602-15





Angled

**For Dilating** 

14.5 cm

Tip Ø: 0.2 mm

No. 18603-14



Curved

Serrated

13 cm

No. 11009-13

15 cm

Tip width: 1.5 mm

No. 11009-15

Straight


1x2 Teeth


12.5 cm

Tip width: 1.3 mm

No. 11024-12







Crile

**Providence-Hospital** 

Straight Straight Curved Curved Curved Straight Curved Serrated Serrated Serrated Serrated Serrated Serrated Serrated 14.5 cm 14 cm 14 cm 14 cm 14 cm 14 cm 14.5 cm Tip width: 1.5 mm Tip width: 1.3 mm Tip width: 1.5 mm No. 13004-14 No. 13005-14 No. 13018-14 No. 13019-14 No. 13013-14 No. 13014-14 No. 13015-14

Adson-Baby

Kelly






**Backhaus Towel Clamp** 












Allis

Tissue and organ holding forceps.



**Lorna Towel Clamp** 





No. 13002-10

No. 13003-10

Curved

**Round tips** 

9 cm


No. 11095-09

### **Clamp Size and Clamp Pressure**

Clamp pressure is determined by the closing force of the clamp (in grams), divided by the area of the vessel wall that is compressed between the jaws (in square millimeters). Thus, the smaller the vessel, the higher the pressure exerted.

Despite the wide range of closing forces, the pressure exerted by each clamp, big or small, remains consistent throughout the series. Each size of clamp exerts a pressure of 5 g/mm<sup>2</sup> when used on the largest vessel in its range, and 15 g/mm<sup>2</sup> when used on the smallest vessel. Even at the higher range, S&T®'s micro clamps have the gentlest working pressure of any small vessel clamp available.

Clamp (ea)	L x W (mm)	Vessel Diameter Range (mm)
00396-01 & 18040-11	3.5 x 1.0	0.4 - 1.0 。 。
00398-02 & 18040-22	5.5 x 1.5	0.6 - 1.5 。 ○
00400-03 & 18040-33	7.5 x 1.75	1.0 - 2.25 。
00325-00	10 x 2.15	2.0 - 5.0









No. 00325-00

Clamp pressure

### Micro Double

For vascular approximation at anastomosis. These clamps are mounted on a 10 mm sliding bar and can easily be moved by applying gentle force. The distance between the jaws of the clamps is adjustable from 3 to 16 mm.



**S&T B-1** 

No. 00396-01



No. 18040-22



No. 18040-33

### **Micro Clamp Applying Forceps**

For use with S&T® vascular clamps.



### Micro Clamps

Compatible with S&T® Clamps and Double Micro Clamps. For clamp sizes B-1, B-2, B-3, and HD-S, as well as Item Nos. 18040-11, 18040-22, and 18040-33.





14 cm Lock: Yes

No. 00071-14

Lock: No No. 00072-14





14.5 cm Lock: No

F.S.T

No. 18040-14



### Micro Clip Applicator with Lock

MOST Popular

# Micro Clip Applicator Forceps Style

Application forceps for Biemer clamps made of Surgical Stainless Steel.

For use with 18058-XX

(Biemer Clip) series

Lock: No


No. 18059-14

**Biemer Clip Applicator Forceps Style** 

### This applicator features special U-shaped side supports for improved seating of the clips.









For use with 18055-XX and 18555-XX series micro serrefines.

Lock: Yes

No. 18056-14





14 cm For use with 18055-XX and 18555-XX series micro serrefines.

Lock: No

No. 18057-14

### **Micro Clips Stainless Steel**

Used to obtain a clear surgical field for extended arterial cross-clamping.



Micro Clips should never be compressed with fingers; use only the appropriate applicator No. 18056-14 or No. 18057-14.

Delicate atraumatic serration



Jaw dimesions (mm)	Total length (mm)	Jaw pressure (g)	Item	
10 x 2	19	60	No. 18055-01	
8 x 2	16	85	No. 18055-02	
6 x 1	15	100	No. 18055-03	
4 x 0.75	13	125	No. 18055-04	

Jaw dimesions (mm)	Total length <sub>(mm)</sub>	Jaw pressure (g)	Item	
6 x 1	17	70	No. 18055-05	
6 x 1	16	70	No. 18055-06	

### **Micro Clips Titanium**

Lightweight with strong clamping pressure.

Delicate atraumatic serration



Jaw dimesions (mm)	Total length (mm)	Jaw pressure (g)	ltem	
10 x 2	19	60	No. 18555-01	
8 x 2	16	85	No. 18555-02	
6 x 1	15	100	No. 18555-03	
4 x 0.75	13	125	No. 18555-04	
	dimesions (mm)  10 x 2  8 x 2  6 x 1	dimesions (mm) length (mm) 10 x 2 19 8 x 2 16 6 x 1 15	dimesions (mm)         length (mm)         pressure (g)           10 x 2         19         60           8 x 2         16         85           6 x 1         15         100	dimesions (mm)         length (mm)         pressure (g)         Item           10 x 2         19         60         No. 18555-01           8 x 2         16         85         No. 18555-02           6 x 1         15         100         No. 18555-03

### **Biemer Clips**

Serrated: Yes Straight

15 mm

18 mm Clamping Length: 6 mm Clamping Length: 9 mm

Closing pressure: 0.29 - 0.39 N / 30 - 40 g



Straight



Biemer clips should never be compressed with fingers; use only the appropriate applicator No. 18059-14.

### Mini Bulldog Artery Clip

Can be used with a minimum of distortion for MRI investigations.



30 mm Serrated: Yes Strong pressure

# Moria

**Moria Vessel Clamp** 

Straight

33 mm Serrated: Yes Light pressure

Jaw dimensions: 13 x 1.5 mm



Moria 201/A

Straight

### Moria De Wecker Clamp

50 mm Serrated: No

Strong pressure Jaw dimensions: 16 x 2.5 mm



### Curved

Moria MC44



No. 18374-44

### Diethrich

Serrated: Yes Closing pressure: 80 g



Straight

No. 18038-45



### Dieffenbach

Serrated: Yes

### Straight 32 mm

Medium pressure



No. 18050-32



Curved

35 mm

Medium pressure

No. 18051-35

Curved

32 mm

Medium pressure

### Straight

35 mm Medium pressure

No. 18050-35

# Straight 50 mm

Strong pressure

No. 18050-50

# Curved

50 mm Strong pressure

No. 18051-50

### **Schwartz**



Serrated: No Strong pressure 26 mm



Straight

No. 18052-01

### Light bend



No. 18052-02



No. 18052-03

### Mini Serrefines

28 mm Medium pressure Serrated: Yes Jaw dimensions: 6 x 2 mm

### Straight



### Curved



No. 18053-28

### Vascular Occluders

### **Applications**

- Fast, dependable zero-flow baseline determination in blood flow studies.
- Partial to full occlusion of vessels for circulation research studies.
- Prolonged implantations, such as chronic blood flow studies.
- Constriction of soft organs in acute or chronic applications.

### **Operation instruction**

- 1. Apply the occluder cuff around the exposed blood vessel. Secure it in place with suture material through the eyelets.
- 2. Occlude the vessel to the desired degree by inflating the diaphragm with air or liquid injected into the actuating tube using a syringe and blunt needle. Clamp the tubing to maintain occlusion of the vessel over time.
- 3. To deflate, simply withdraw the air or liquid.

### **Benefits**


- Occlusion is accomplished without traction on vessels or surrounding
- The actuating tube may be exteriorized from the occlusion site for remote actuation.
- Performance during implantations is reliable and consistent.
- Easily maintained: can be autoclaved or cold sterilized.
- Made from soft, flexible 100% silicone rubber.
- Fully operational with instant response using air, liquids, or inert gases.

### Selection

The lumen diameter when deflated is usually the determining factor in the selection of the proper device. Select the occluder size that provides a slightly loose fit around the subject vessel. The cuff width and thickness will vary, depending on models and sizes.



The actuating tube length is 90 cm and measures 0.76 mm ID x 1.6 mm OD. Detailed instructions are included.



A: Lumen diameter

B: Cuff width

C: Cuff thickness

VO-1.5 N A: 1.5 mm - B: < 3.5 mm - C: 1.5 mm

No. 18080-01

V0-2

A: 2 mm - B: 5 mm - C: 1.5 mm No. 18080-02

VO-3 A: 3 mm - B: 5 mm - C: 1.5 mm

No. 18080-03

### V-Hook

Aids in feeding fine suture around blood vessels, works very similar to ligation aids.

Magnified View

### **Ligation Aids**


A handy instrument for leading a fine suture under and around a blood vessel for tying-off purposes. Also useful for ligating peripheral nerves. Grooves toward hole help threading.



Magnified View















### **Vessel Dilators**

Highly modified forceps with elongated, parallel tips which are highly polished. Designed for controlled intraluminal vessel dilation.

S&T®







3)

11 cm Tip Ø: 0.1 mm

No. 00126-11







1

Angled 45°

13.5 cm Tip Ø: 0.2 mm No. 00276-13



Angled 10°

12 cm

Tip Ø: 0.1 mm No. 18131-12

Tip Ø: 0.2 mm No. 18132-12


This micro cannula is most helpful for quick and easy insertion in arteries and veins. It can be used as an "indwelling" catheter during long term experiments. With care, thorough rinsing and sterilization, it can be reused several times. Detailed instructions are included.

### Features

- Biocompatible Teflon® tubing.
- Cannula size 0.33 mm OD 0.16 mm ID
- A pointed Surgical Stainless Steel stylus (trocar), to aid in vascular penetration and allow a maximum range of cannula manipulation during insertion.
- Contoured tip for easy insertion, fitted with shoulder for "tie-in" to prevent accidental removal.
- Female luer connection for linkage to infusion lines.



Cannulation Forceps (Item No. 18403-11) and Spinal Cord Hook (Item No. 10162-12) aid in the insertion process.



### 2 Insertion

Remove trocar and advance cannula.



### 3 Secure

Use shoulder hub and suture to "tie in" cannula to prevent accidental removal.



Instruments sold separately.

8 cm Cannula size: 0.33 mm OD - 0.16 mm ID

No. 18000-10

### **Tubing Introducer Forceps**

These vessel cannulation forceps are useful whenever a fine plastic tube has to be introduced into a small blood vessel of almost equal size. The hollowed jaws hold the tubing securely without deforming it. The tip of the tubing is directed exactly into the vessel opening without damage to the vessel from unwanted movement.







For tubing of 0.35 mm OD

11 cm

No. 00608-11


For tubing of 0.6 mm OD

11 cm

No. 00574-11







No. 11280-11

For tubing of 0.5 mm OD 13 cm

No. 11281-11

For tubing of 1 mm OD 13 cm

No. 11282-11

For tubing of 0.35 mm OD 11.5 cm

No. 18403-11

For tubing of 0.6 mm OD 11.5 cm

No. 18406-11

For tubing of 0.6 mm OD 13.5 cm

No. 18406-13

# PROBES, NEEDLE HOLDERS, WOUND CLOSURE, RETRACTORS


Probes & Hooks	133 – 138
Spatulae	139 — 142
Spoons	143
Pins & Holders	144 – 145
Wound Closure	146 — 156
Retractors	157 — 163

### **Double Pronged Hook**

For peritoneum, pleura, endocardium, organ capsules and skin. Grasps lifts skin, tissue and dura etc.



Magnified View




The tip consists of two very sharp, nearly horizontally opposed points. By applying gentle downward pressure and twisting slightly in a clockwise fashion, the points will embed themselves horizontally into fine membranes or skin etc., so that it may be lifted upward. This instrument also has been used successfully when separating fine tissue from an underlying area that is not to be cut.



**Spinal Cord Hook** 

Magnified View




Shaped like a ligation aid and sized to isolate the spinal cord. The cord can then be cut with minimal damage to the dura.



No. 18067-11



No. 10162-12



These instruments are lightweight, color coded for easy identification, have an ergonomic design and a non-reflective surface.



Micro Point



12.5 cm Tip Ø: 0.4 mm

No. 10165-11



Y-Hook







·TITANIUM ·

Magnified View

Micro Spatula


Angled 30°

12.5 cm

Tip: 0.25 x 0.6 mm

No. 10167-11





**Gross Anatomy Probe** 





**Grooved Director** 

The Grooved Director is a specialized multipurpose device which is used by surgeons in different areas to probe cavities and wounds for example.



13 cm with Tongue Tie

No. 10290-13

F·S·T

PROBES & HOOKS

Angled

12.5 cm

Tip Ø: 0.6 mm


No. 10160-13

Angled 45°

15 cm

Tip Ø: 1.2 mm


No. 10088-15



Kit with 4 pieces.












**PRICES SUBJECT TO CHANGE** PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER









No. 10098-06

**Mirror Finish Cutting** 

Spoon

Non stick cutting spoon.

### **Mini Hippocampal Dissecting Tool**

The instrument ends are malleable to allow you to adjust their shape to the correct angle for optimal positioning.

### **Hippocampal Dissecting Tool**

The instrument ends are malleable to allow you to adjust their shape to the correct angle for optimal positioning.

### Spoon

Moria

Moria 1121B

12.5 cm

Tip Ø: 8 mm


Depth: 2 mm

No. 10321-08

**Mini Perforated Spoon** 

Most commonly used to take embryos out of solution and handling very small organs.

### Moria



Most commonly used to take embryos out of solution and handling very small organs.
















Tip thickness: 0.5 mm No. 10099-15








10 cm Tip Ø: 8 mm Depth: 1 mm Hole Ø: 0.8 mm






Moria MC 17 BIS



Moria MC 17

14.5 cm Tip Ø: 20 mm Depth: 5 mm Hole Ø: 0.9 mm

No. 10360-13 No. 10370-17



Tip thickness: 0.5 mm

No. 10099-12



PRICES SUBJECT TO CHANGE PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER

13 cm

Tip Ø: 4.5 mm

Depth: 0.15 mm

### Pin/Needle Holder

Light weight.






17 cm Aluminum handle with nickel-plated chuck Jaw opening: 0 to 1 mm

No. 26018-17

### Moria Pin/Needle Holder







Moria MC1


12 cm Nickel-plated Jaw opening: 0 to 1 mm

No. 26016-12

### **Needle Inserts**

All insert pins fit any of our pin holders. Pins are antimagnetic and acid resistant.

Surgical Stainless Steel (Inox)



### **Minutien Pins Black Anodized**

These high quality pins are black for best visibility under high magnification and bright light.

> Tip Ø: 0.0175 mm Rod Ø:0.15 mm 1.2 cm

Box of 500 pins No. 26002-15

### **Minutien Pins Stainless Steel**

The thinnest and finest pins available.



Tip Ø: 0.0125 mm Rod Ø: 0.1 mm 1.2 cm

Box of 500 pins No. 26002-10 Tip Ø: 0.0175 mm Rod Ø: 0.15 mm 1.2 cm

Box of 500 pins No. 26003-15 Tip Ø: 0.02 mm Rod Ø: 0.2 mm 1.2 cm

Box of 500 pins

No. 26002-20

### **Insect Pins Black Anodized** 5.5 cm Nylon heads Tip Ø: 0.03 mm 4 cm Pkg. of 100 pins 2 Size 00 0 5 1 3 4

-25 -30 -35 -40 -45 -50 -55 -60 -65 -70

### **Insect Pins Stainless Steel**

Nylon heads Tip Ø: 0.03 mm

Rod Ø (mm)

No. 26000

Pkg. of 100 pins

Size 0 1 2 3 4 5 6 Rod Ø (mm) No. 26001 -25 -30 -35 -40 -45 -50 -55 -60 -65 -70

### **Etched Tungsten Micro Needles**

These Tungsten micro needles are drawn to an approximate 1 micron tip and three different rod diameters are offered. Each provides lateral stability as well as minimal spring action during manipulation. The micro needles can be used with our pin holders shown on page 144.





Tip Ø: 1 μm Box of 10 needles

(mm)	item
0.125	No. 10130-05
0.25	No. 10130-10
0.5	No. 10130-20

F·S·T

PINS & HOLDERS

Micro pattern, with slightly blunted tips.







Never use a micro needle holder with needles that are too thick, as this can easily damage the fine jaws. Use only micro needles of less than 0.2 mm diameter (see page 152).



**WOUND CLOSURE** 

9 cm Tip width: 0.6 mm

Lock: Yes

No. 12060-01

Tip width: 0.6 mm Lock: No

No. 12060-02



9 cm Tip width: 0.6 mm Lock: Yes

No. 12061-01

Tip width: 0.6 mm Lock: No

No. 12061-02

### Round-Handled Micro/Macro

These Round-Handled Needle Holders have a dual-purpose and can be used for microscopic and macroscopic applications. The fine tapered tips can be used to handle small needles up to 0.2 mm in diameter. While the remaining portion of the tip ensures a secure hold for larger needles. The balanced design improves control and stability, while the ergonomic handle enhances grip and reduces hand strain.







Tip width: 0.4 mm Lock: No

No. 12875-13



Round handles allow easy fingertip adjustments.

### Round Handled with Suture Cutter

Round handles allow easy fingertip adjustments.





12.5 cm Lock: Yes

No. 12075-12 No. 12074-12

Lock: No

Lock: No

No. 12076-12



Curved

14 cm Cutting edge: 3.5 mm Lock: No

No. 12075-14

### Kalt

The Kalt needle holder is characterized by its extraordinary design and good handling. The closure system closes softly and securely and can be easily opened again with your thumb.







traight

Serrated: Yes 14 cm Lock: Yes

No. 12080-14



Curved

Serrated: Yes 14 cm Lock: Yes

No. 12081-14

### Castroviejo with Tungsten **Carbide Jaws**







Straight

14 cm Tungsten Carbide jaws Lock: Yes

No. 12565-14

### Mathieu

Extra slim and narrow profile, with squeeze-lock. Facilitates excellent wrist motion for suturing.







14 cm





Tungsten Carbide jaws Lock: Yes

No. 12010-14




14 cm

No. 12510-14









13 cm Smooth Lock: Yes No. 12001-13



13 cm Tungsten Carbide jaws Lock: Yes

No. 12501-13

### Halsey Micro with **Tungsten Carbide Jaws**











12 cm Tungsten Carbide jaws Lock: Yes

No. 12500-12

F·S·T

**WOUND CLOSURE** 





LEFT-HANDED

TUNGSTEN

CARBIDE

12502-12

12512-12



Crile-Wood

### Slim Elongated with Tungsten **Carbide Jaws**





12504-15



TUNGSTEN CARBIDE





12 cm Cutting edge: 10 mm Lock: Yes

No. 12002-12

No. 12102-12



14 cm

Cutting edge: 15 mm Lock: Yes


No. 12002-14





Tungsten Carbide jaws Cutting edge: 10 mm Lock: Yes No. 12502-12

No. 12512-12



Tungsten Carbide jaws Cutting edge: 15 mm Lock: Yes No. 12502-14



Lock: Yes

No. 12004-16





No. 12004-18

70% of

actual size



15 cm

Lock: Yes

No. 12003-15







15 cm

Lock: Yes No. 12504-15



15 cm Tungsten Carbide jaws

Lock: Yes

No. 12005-15

### **Microsurgical Needles with Suture Thread Attached**


The combination of proprietary needle design technology and state-of-the-art manufacturing techniques has enabled FST® to provide a microsurgical needle that is unmatched in the industry. FST microsurgical needles provide frictionless tissue penetration. They are inspected under the microscope at critical stages during the manufacturing process, and a multi-point final inspection of each suture and needle combination ensures a high quality product. The suture-to-needle attachment is tested to higher standards than the USP specifications.

### Features:

- Needle and thread combination comprising of nonabsorbable black nylon polyamide monofilament thread and Surgical Stainless Steel 3/8 circle needle. Nylon is suitable for closing ophthalmic, neurological and microvascular wounds.
- 3/8 circle precision taper point. Surgical Stainless Steel needle honed under microscope by hand. Superfine round tip with flattened body to improve stability in needle holder. Easy penetration. Precise stitch placement.
- Grinding geometry and suture needle attachment are directed toward frictionless tissue penetration.
- Each needle suture combination is individually packed and sealed in a peel-open blister pack.

Sterile

Pkg. of 12 disposable units






400% of actual size



No. 12052-10





Suture size	Suture length (cm)	Needle Ø (micron)	Needle Ø	Needle chord length (mm)	Circle	Item
8/0	15	140	0.14	29	3/8	No. 12052-08
9/0	15	100	0.1	24	3/8	No. 12052-09
10/0	10	70	0.07	18	3/8	No. 12052-10
11/0	10	50	0.05	16	3/8	No. 12052-11


### **Kalt Suture Needles**



- Polished Surgical Stainless Steel
- 3/8 curved
- Round with standard eye
- Suitable for sutures up to 0.32 mm diameter Ø: 0.35 mm

Sterile: No

Pkg. of 12 needles



MOST Popular

Size	Needle Chord Length	Item
3	12 mm	No. 12050-03
2	16 mm	No. 12050-02
1	20 mm	No. 12050-01
0	23 mm	No. 12050-10
00	26 mm	No. 12050-20

### **Silk Suture**

Silk suture has a smooth flow through tissue while maintaining knot security.

Braided Silk

### Sterile: No

Can be ETO sterilized



Size	Suture Ø (mm)	Spool length	Item
7/0	0.05	91	No. 18020-70
6/0	0.08	23	No. 18020-60
5/0	0.12	23	No. 18020-50
4/0	0.17	23	No. 18020-40
3/0	0.22	23	No. 18020-30
2/0	0.32	23	No. 18020-20
0	0.37	23	No. 18020-00
1	0.45	23	No. 18020-01
2	0.55	23	No. 18020-02
3	0.65	23	No. 18020-03

### For Suture Tying Forceps, see pages 99 - 101 und 111.

### **Nylon Monofilament Suture**

Nylon is known for its uniform diameter flexibility and outstanding consistent knot security.

Nylon monofilament

### Sterile: No

Can be ETO sterilized



Size	Suture Ø	Spool length	Item
8/0	0.04	91	No. 18030-80
7/0	0.05	91	No. 18030-70
6/0	0.07	91	No. 18030-60
5/0	0.1	91	No. 18030-50
4/0	0.15	91	No. 18030-40
3/0	0.2	91	No. 18030-30

# **WOUND CLOSURE**

### **Michel Suture Clips**

Anti-magnetic

WOUND CLOSURE

- Acid-resistant Surgical Stainless Steel
- Pkg. of 1,000 clips

Size (mm)	Item
7.5 x 1.75	No. 12040-01
11 x 2	No. 12040-02
12 x 2.5	No. 12040-03



# **Applying Forceps** F-S-T

12.5 cm For clip sizes 12040-02 and 12040-03


No. 12018-12

## **Applying and Removing Instruments** Angled Curved 13 cm 13 cm For clip size 12040-01 For clip sizes 12040-02 and 12040-03 No. 12029-12 No. 12028-12

### **Reloadable Skin Stapler with Rotating Head**

The 360° swivel nose and unique design offers maximum visibility at the delivery point. The Surgical Stainless Steel staples are non-sterile and feature a patented staple design which make for quick and easy removal. The disposable staple cartridge includes a printed staple counter.

### The handle is reusable and autoclavable. Staples are non-sterile. Components sold separately or as a kit.



Compatible with Reloadable Skin Stapler with Rotating Head and Manipler® Skin Stapler.

**Surgical Staple Remover** 



10 cm

No. 12030-04

### Manipler® Skin Stapler

Quality, slim, lightweight design with consistent actuation and staple count indicator.




Fomon

### AutoClip® System

The AutoClip® System is a safe, secure method of wound closure. Faster and easier than suture, this system is perfect when staples may be too thin to provide adequate closure. Non-sterile. Components sold separately or as a kit. We recommend using wet steam sterilization at 120-125 °C for at least 30 minutes in an autoclave to sterilize the device.







Clips



Remover


13 cm 2 Blunt

No. 17021-13

Guthrie


Extra delicate

pattern.



No. 17023-13

Fine




F-5-T 17022-13

Blair







2 with 2 mm Ball each tip

No. 17020-17

### AutoClip® Kit

1 Applier

WOUND CLOSURE

- 9 mm clips (100)
- 1 Remover
- Storage case

No. 12020-00

AutoClip® applier

**Applier** 

11.5 cm

No. 12020-09

9 mm Pkg. of 100 clips

No. 12022-09

AutoClip® clips

9 mm

Pkg. of 1,000 clips No. 12022-19

AutoClip® clips


No. 12023-00

AutoClip® remover

12.5 cm

F·S·T

### **Goldstein Mini with Blades**







2 cm Max. spread: 1 cm Atraumatic blades

No. 17004-03

RETRACTORS

### Goldstein Mini





2 cm Max. spread: 0.5 cm 3 x 3 Sharp teeth

No. 17002-02

### **Barraquer Mini**







2 cm Max. spread: 1 cm

No. 17000-01

### **Corcelle Speculum**



Moria 7699

4.5 cm Max. spread: 3.5 cm Curved

No. 17370-53

### **Fine Adjustable Spring Wire**




5 cm Max. spread: 3.5 cm Curved

No. 17004-05

### **Goldstein Mini with Blunt Teeth**






2 cm Max. spread: 1 cm 2 x 2 Blunt teeth

No. 17005-03

### Goldstein






3 cm Max. spread: 1.8 cm 3 x 3 Sharp teeth No. 17003-03

### Barraquer





3 cm Max. spread: 1.5 cm

No. 17000-03

4 cm



3 cm Max. spread: 1.2 cm

No. 17000-04

No. 17000-02

### Barraquer







Moria 6346D

3.4 cm Max. spread: 1.5 cm

No. 17300-02



Moria 6346C

3.3 cm Max. spread: 1.5 cm

No. 17300-03



Moria 6346B

3.5 cm Max. spread: 2 cm

No. 17300-04

### Agricola





4 cm Max. spread: 3.5 cm 3 x 3 Sharp teeth

No. 17005-04









7 cm Max. spread: 4 cm WireSlightly curved down

No. 17010-07



7 cm Max. spread: 4 cm Slightly curved 4 x 4 Blunt



7 cm Max. spread: 4 cm Blunt teeth 4 x 4 Blunt

No. 17009-07



8 cm Max. spread: 4.5 cm Blades

No. 17009-08



6 cm Max. spread: 3 cm

14 x 10 mm

No. 17006-05



Max. spread: 5 cm

22 x 16 mm

No. 17007-08

No. 17008-07

### Heiss

### Weitlaner







4 x 4 Blunt teeth

10 cm Max. spread: 2.3 cm

No. 17011-10



Straight 2 x 3 Blunt teeth

Straight 3 x 4 Blunt teeth

10.5 cm Max. spread: 4 cm

13 cm Max. spread: 4 cm

No. 17012-11

No. 17012-13



### **Magnetic Fixator Retraction System**

This Small Animal Retraction System simplifies set-up procedure and provides the examiner with an excellent overview of the operative site. It consists of three components: Base Plate (Nos. 18200-03 & 18200-04); Fixators (Nos. 18200-01 & 18200-02); and Retractors (Nos. 18200-05 through 18200-12).

The Base Plates are made from magnetic Stainless Steel and are available in two sizes. The surface is nonreflective and both plates have a window in the center which allows the animal contact with body temperature control systems. The window sizes are 10 x 13 cm for the small base plate and 15 x 25 cm for the large.

The Fixators are movable anchors made from anodized aluminum with a surgical grade rare earth magnet in the base for attachment anywhere on the Base Plate. The magnets do not interfere with lab electronics and can be used in close proximity (3 cm) to ferro-magnetic items. They incorporate pushbutton-activated spring locks that grip either the Wire or Elastomer Retractors.

The Retractors are interchangeable Retractor Tips (Nos. 18200-08 through 18200-12) which attach to either Retractor Wire (Nos. 18200-05 & 18200-06) for static retraction or to Elastomer (No. 18200-07) for dynamic retraction. The Retractor Tips are light and atraumatic to delicate tissue or membrane. All items are autoclavable, except for the elastomers.

### Kit for animals up to 200 g

Includes:

- 1 Small base plate (No. 18200-03)
- 2 Short fixators (No. 18200-01)
- 1 Tall fixator (No. 18200-02)
- 1 Elastomer, 2 m roll (No. 18200-07)
- 1 Package of retractors, blunt, 1 mm wide (No. 18200-09)
- 1 Package of retractors, blunt, 2.5 mm wide (No. 18200-10)

No substitutions No. 18200-20



### Kit for animals from 200 g up to 5 kg

Includes:

- 1 Large base plate (No. 18200-04)
- 2 Short fixators (No. 18200-01)
- 2 Tall fixators (No. 18200-02)
- 1 Elastomer, 2 m roll (No. 18200-07)
- 1 Package of retractors, blunt, 2.5 mm wide (No. 18200-10)
- 1 Package of retractors, blunt, 5 mm wide (No. 18200-11)

No substitutions

No. 18200-50

### Individual parts:

Short fixator, 3 cm high

No. 18200-01

Tall fixator, 6 cm high

No. 18200-02

Small base plate, 20 x 30 cm Each

No. 18200-03

Large base plate, 25 x 35 cm Each

No. 18200-04

Wire handles, 10 cm long Pkg. of 10

No. 18200-05

Wire handles, 14 cm long Pkg. of 10 No. 18200-06

Elastomer, 2 m roll

No. 18200-07

Sharp retractors, 1 mm wide Pkg. of 10

No. 18200-08

Blunt retractors, 1 mm wide Pkg. of 10

No. 18200-09

Blunt retractors, 2.5 mm wide Pkg. of 10

No. 18200-10

Blunt retractors, 5 mm wide Pkg. of 10

No. 18200-11

Blunt retractors, 7.5 mm wide Pkg. of 10

No. 18200-12

### ACCESSORIES & INSTRUMENT CARE

Animal Accessories	165 — 170
Instrument Care & Accessories	171 – 182
Surgical Accessories	183 – 186
Laboratory Accessories	187 – 197



### Reusable Feeding Needles Ball Tip

Surgical Stainless Steel feeding needles with luer lock hub. Also ideal for transcardiac perfusions with a fixative. The needles come with a round tip.



### No. 18061-24

Gauge	Length (mm)	Tip Ø	Shape	Mice (g)	Rats (g)	Guinea pigs	Item
24	25	1.25	Straight	>14			No. 18061-24
22	25	1.25	Straight	15-20			No. 18061-22
20	38	2.25	Straight	20-25	50-75		No. 18064-20
20	38	2.25	Curved	20-25	50-75		No. 18065-20
18	38	2.25	Straight	20-30	100-200	250-300	No. 18061-38
18	50	2.25	Curved	25-30	100-200	250-300	No. 18061-50
16	75	3	Curved		150-300	350-450	No. 18061-75
16	100	3	Curved		150-300	350-450	No. 18061-10

### These needles come with a barrel tip.






### No. 18060-20

Gau	ge	Length (mm)	Tip Ø	Shape	Mice (g)	Rats (g)	Guinea pigs	Item
	24	25	1.7	Straight	>14			No. 18060-24
	22	25	1.7	Straight	15-20			No. 18060-22
	20	30	1.7	Straight	20-25	50-75		No. 18060-20
	20	38	1.7	Straight	20-25	50-75		No. 18063-20
	20	38	1.7	Curved	20-25	50-75		No. 18067-20
	18	38	1.7	Straight	20-30	100-200	250-300	No. 18060-18
	18	50	2	Curved	25-30	100-200	250-300	No. 18067-50
	16	75	2.4	Curved		150-300	350-450	No. 18067-75
	16	100	2.4	Curved		150-300	350-450	No. 18067-10

### **Single Use Feeding Needles**

These economical feeding needles feature plastic luer hubs with silicone tips. They are nontoxic and nonpyrogenic. Needles can be easily formed to any desired shape. Sterile.






G	auge	Length (mm)	Tip Ø	Mice (g)	Rats (g)	Guinea pigs	Item
	20	38	1.75	20-25	50-75		No. 18061-20
	18	50	2.5	25-30	100-200		No. 18061-18
	15	75	3		200-350	400-600	No. 18061-15

### **Disposable Flexible PTFE Feeding Needles**

These affordable needles feature flexible PTFE tube with silicone tip and plastic luer hub, delivered in sterile pouches. They are very flexible, and translucent for easy confirmation of reagent delivery. The soft silicone tips are non-toxic and non-pyrogenic and ensure that no damage occurs during the insertion. Sterile.

### Pkg. of 5 needles






Gauge	Length (mm)	Tip Ø	Mice (g)	Rats (g)	Guinea pigs	ltem
20	38	1.8	20-25	50-75		No. 18961-20
18	38	2	20-30	100-200	250-300	No. 18961-18
17	52	2.4	25-30	100-200	250-300	No. 18961-17
15	78	2.6		200-350	400-600	No. 18961-15
18	75	2	25-30	150-250	250-300	No. 18961-75
17	85	2.4	25-30	150-250	250-300	No. 18961-85

### GenoPunch

The interchangeable GenoPunch aids in the fast and simple collection of the 2 mm tissue sample for DNA analysis and/or genotyping. FST® Exclusive Item. Autoclavable.







- 1 GenoPunch
- 1 0.5 ml tube holder
- 1 1.6 ml tube holder
- 1 Screwdriver
- 4 Screws

FST does not supply microtubes.



11.5 cm Surgical Stainless Steel (Inox) and anodized aluminum Punch Ø: 2 mm

No. 24215-02

The Ear Puncher relieves stress on the wrist during ID marking as it requires only a very gentle pressure to cut a hole. The Stainless Steel puncher has a sloped edge that provides a cut rather than a punch, which minimizes discomfort. The ear puncher should be lubricated to improve longevity. Autoclavable.



14 cm Punch Ø: 1 mm

14 cm Punch Ø: 2 mm

No. 24916-01

No. 24916-02

F·S·T

PRICES SUBJECT TO CHANGE PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER

Designed to ensure a fast and clean punch, this ear punch is lightweight and fits comfortably in the palm of your hand.

Autoclavable.





11.5 cm Surgical Stainless Steel (Inox) and anodized aluminum Punch Ø: 2 mm

No. 24214-02

### **Finger Loop Ear Punch**

This ear punch provides stability and control for researchers using the numbering system. The 1 mm ear punch offers greater numbering variations and/or patterns.



10.5 cm Punch Ø: 2 mm

No. 24212-02

Punch Ø: 1 mm

No. 24212-01

### Marking Pen

The Marking Pen was developed for easy, long lasting identification of laboratory animals. The markers contain no harmful chemicals and instead nontoxic pigment approved by the Japanese Ministry of Health for use in food, drugs and cosmetics. It is intended for use on the fur rather than the skin and will last for 6-12 weeks.



Color can stain skin and clothing. The use of gloves and lab coats is recommended. These colors will not work on black fur.



Green Fine Tip 4.5 x 7.5 mm 10 ml No. 24202-01

No. 24202-02

Fine Tip Fine Tip 4.5 x 7.5 mm 4.5 x 7.5 mm 10 ml 10 ml No. 24202-03

Yellow

Blue

Fine Tip 4.5 x 7.5 mm 10 ml No. 24202-04

Red

Fine Tip 4.5 x 7.5 mm 10 ml

No. 24202-05

Orange

No. 24202-15

Orange

Fine Tip 4.5 x 7.5 mm 10 ml

Black

No. 24202-06



Green Standard Tip 15 x 10 x 10 mm

15 ml No. 24202-11

Standard Tip 15 x 10 x 10 mm 15 ml

No. 24202-12

Blue

Standard Tip 15 x 10 x 10 mm 15 ml

No. 24202-13

Yellow

Standard Tip 15 x 10 x 10 mm 15 ml

Red

No. 24202-14

Standard Tip 15 x 10 x 10 mm 15 ml

Standard Tip 15 x 10 x 10 mm 15 ml

Black

No. 24202-16

### Skin Marker

This marker features a versatile tip for precise fine (0.2 mm) line marking, or when held at an angle, a broad (2 mm) line marking. The marker's violet ink is nontoxic, non-tattooing and surgical prep-resistant.

Markings can be removed from the skin with isopropyl alcohol. A six-inch ruler calibrated in both centimeters and inches is included.



No. 18000-30

50% of actual size

### **Self Piercing Ear Tags**

Ideal for use in small animals such as mice and rats, these ear tags come numbered from 1-100 or 1-500. Silver colored tags with black numbers provide high visibility making them easy to read. Once applied, tags seal and cannot be removed, providing an easy to use, permanent identification method. These ear tags are safe for use in salt water.






### Custom alphanumeric sequences and numbers above 1000 are available by special order.

Customized ear tags available in packages of 100 or 500. The following markings are possible: 1 to 6 numbers only or 1 to 6 letters and numbers.

For special order, use item No. 70000-10 and indicate sequence. Please allow 4 to 5 weeks for delivery.




### **Fine Forceps Repair Kit**

Realign and repair fine forceps with this convenient kit. Bent forceps tips can be re-straightened with the chain nose pliers. The sharpening stone has two grooves, 1.6 mm and 0.8 mm, to assist in resharpening the forceps tips. O-rings clamp the forceps shut to ensure that both shanks are equally aligned and sharpened. Instrument oil assists in producing a smooth sharp edge. Detailed instructions are included.



### Kit includes:









### **O-Rings for Forceps**

These rubber O-rings fit our forceps and keep them closed with variable degrees of pressure on the tips.

Pkg. of 10

No. 11200-00





### **Sharpening Stone for Fine Forceps**

This natural Arkansas stone assists in resharpening the tips of fine forceps. The stone has two grooves, 0.8 mm and 1.6 mm and comes in a cedar box with ten O-rings.

12.5 x 6.5 x 3 cm

F·S·T

50 ml bottle No. 29055-00

Oil bottle or spray

Physiologically inert oil for lubricating

surgical instruments, box locks, ratchets

and any overlapping joints.

No. 29008-22

Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 8 colors available. **Autoclavable.** 



Color	Roll length	Tape width	Material	Item
Green	290	0.6	Vinyl	No. 29071-01
Blue	290	0.6	Vinyl	No. 29071-02
Yellow	290	0.6	Vinyl	No. 29071-03
Red	290	0.6	Vinyl	No. 29071-04
White	290	0.6	Vinyl	No. 29071-05
Black	290	0.6	Vinyl	No. 29071-06
Orange	290	0.6	Vinyl	No. 29071-08
Purple	290	0.6	Vinyl	No. 29071-09

0.6

No. 29080-02

### Instrument Marking Tape 7 assorted colors

Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. Includes 7 assorted color rolls: Yellow, Red, Black, White, Green, Blue and Orange. **Autoclavable.** 

107





### Instrument Marking Tape - 510 cm

Organize instruments quickly and easily with this economical, non-toxic, pressure sensitive tape. Color is unaffected by boiling, cold sterilization or autoclaving. 9 colors available. **Autoclavable.** 



Color	Roll length	Tape width	Material	Item
Green	510	0.6	Vinyl	No. 29070-01
Blue	510	0.6	Vinyl	No. 29070-02
Yellow	510	0.6	Vinyl	No. 29070-03
Red	510	0.6	Vinyl	No. 29070-04
White	510	0.6	Vinyl	No. 29070-05
Black	510	0.6	Vinyl	No. 29070-06
Brown	510	0.6	Vinyl	No. 29070-07
Orange	510	0.6	Vinyl	No. 29070-08
Purple	510	0.6	Vinyl	No. 29070-09

### **Self Sealing Sterilization Pouches**

Preformed pouches with self adhering tape are blue tinted on one side for visualization of contents and paper-backed on the other side. Pouches will work with either steam autoclave or ethylene oxide sterilization.

### Material properties:

- Protects instruments during sterilization.
- Provides sterile storage until use.
- Write on pouches for easier identification during sterilization.



Outer: 9 x 17 cm | Inner: 7 x 11.5 cm For instruments up to 11 cm long

Pkg. of 100 pouches

No. 28000-06

Outer: 13.5 x 28 cm| Inner: 11.5 x 22.5 cm For instruments up to 20 cm long

Pkg. of 100 pouches

No. 28000-11

Outer: 19 x 36 cm | Inner: 17 x 30 cm For instruments up to 28 cm long

Pkg. of 100 pouches

No. 28000-14

The sterilizer provides a safe, effective and convenient method for sterilizing your surgical instruments prior to and during surgery or tissue culture. Without using potentially dangerous gases flames or chemicals, you can ensure that your instruments will be free of pathogens and microbial contaminants.

The beads take 15–20 minutes to reach 240 °C to 270 °C (464 °F to 518 °F), then the sterilizer is ready for use. The temperature is thermostatically maintained and the unit can be left on all day. Instruments require approximately 20 seconds of contact with the beads for complete sterilization. The sterilizer is safe to use even in laminar flow hoods.

Sterilization test results conducted by the Pasteur Institute in Lille, France, are available upon request. Available in 110 V and 220 V. CE Certified and UL Approved.

### Sterilizer includes:

- Sterilizer unit
- 1 Bag of glass beads
- 1 Dust cover for well
- 1 Small heat-tempered glass dish for small objects



250 shown.

Instruments not included.

60% of actual size

### 250 Sterilizer

Well Ø: 4 cm Well depth: 8 cm Well volume: 100 ml Total unit size: 14.5 x 12.5 x 13.5

No. 18000-45

### Replacement glass beads for 18000-45

Weight: 150 g

No. 18000-46

### 350 Sterilizer

Well Ø: 4 cm Well depth: 14 cm Well volume: 175 ml Total unit size: 20.5 x 12.5 x 13.5

No. 18000-50

### Replacement glass beads for 18000-50

Weight: 300 g

No. 18000-51

### **Micro Bead Sterilizer**

Micro Bead Sterilizer is ideal for the sterilization of small and large research tools (forceps, scissors, etc.). When the chamber is filled with the included glass beads, high temperatures (up to 300 °C / 572 °F) can be used to eliminate the presence of bacteria, spores and other micro-organisms.

### Small Sterilizer 115 V

Well Ø: 4 cm Well depth: 8 cm Well volume: 100 ml Total unit size: 14.5 x 15.5 cm

No. 18090-45

### Small Sterilizer 220 V

Well Ø: 4 cm Well depth: 8 cm Well volume: 100 ml Total unit size: 14.5 x 15.5 cm

No. 18090-46

### Replacement glass beads for small sterilizer

Weight: 125 g

No. 18000-47



### Tall Sterilizer 220 V

Well Ø: 4 cm Well depth: 14 cm Well volume: 175 ml Total unit size: 14.5 x 15.5 x 21.5 cm

No. 18090-51

### Replacement glass beads for tall sterilizer

Weight: 256 g

### Tall Sterilizer 115 V

Well Ø: 4 cm Well depth: 14 cm Well volume: 175 ml Total unit size: 14.5 x 15.5 x 21.5 cm

No. 18090-50


### **Autoclave Tape**

Tape will permanently indicate diagonal stripes after exposure to steam autoclave temperature between 121 °C and 130 °C during a 40 - 60 minute time period. Use with sterilization pouches, sterilization wraps and sterilization cases.

### The tape can be used in:

- Gravity sterilizers operating at 121 °C (+3 °C) for 30 minutes or
- Pre-vacuum sterilizers operating at 132 °C (+3 °C) for 4 minutes and 135 °C (+3 °C) for 3 minutes
- Shelf life/durability is 36 months





1,90 cm x 55 m

No. 28001-02

### **Disposable Sterilization Wraps**

These biodegradable cellulose fiber wraps are resistant to tears and conform easily to maintain a breathable yet impervious barrier to airborne bacteria, saline, blood and water. Wrapping your instrument case securely will ensure the sterility of your instruments until ready for use. Available in two sizes.



Cases not included.

38 x 38 cm Pkg. of 25 wraps

No. 28015-15

51 x 51 cm Pkg. of 25 wraps

No. 28020-20

### **Silicone Mat**

High quality, extra soft mat for safe and secure placement of delicate instruments. This mat can easily be cut along the premarked grooves to fit in various sized instrument cases. Autoclavable.



30 x 20 cm

No. 20325-00

### **Magnetic Sterilization Mat**

Used during sterilization and for safe anti-slip placement around surgical settings. Autoclavable.




Instruments not included.

30 x 20 cm

No. 20350-00

### Silicone Mat

Large, economical silicone mat. Autoclavable.




46 x 25 cm

No. 20325-10

# INSTRUMENT CARE & ACCESSORIES

### **Plastic Sterilization Container**

This plastic sterilization container includes a silicone mat for safe placement of instruments. Perforated top is clear for visual identification of instruments. Suitable for autoclave or ethylene oxide sterilization.



### **Stainless Steel Sterilization Container**

This container will assist in the sterilization and storage of your fine surgical instruments. It comes with a custom fit silicone mat to keep instruments in place while the sterilizing process and storage. The textile filters are necessary for keeping instruments sterile while they remain in the container. Suitable for autoclave or ethylene oxide sterilization.



30 x 14 x 4 cm No. 20850-00 5 Long-term textile filters No. 20850-10


1 Long-term textile filter

No. 20850-11

### **Plastic Sterilization Container**

This sterilization container is made from Ultem® polymer to ensure strength and long lasting durability. The transparent perforated lid provides easy identification of your surgical instruments. Silicone mat is included for a safe and secure hold. Trays are also stackable for convenience. Suitable for autoclave or ethylene oxide sterilization.








19 x 7 x 2 cm Single tier

No. 20840-01

17 x 10 x 2 cm Single tier

No. 20840-02





25 x 15 x 4 cm Double tier

No. 20842-03

25 x 15 x 2 cm Single tier

No. 20840-03

### **Colored Aluminum Instrument Cases**

This instrument case comes with a spiked silicone mat for secure storage of surgical instruments. Five colors available but may vary due to lot to lot variations.

The case and silicone insert are autoclavable but should not be used for sterilization of surgical instruments.



### **Moria Instrument Case**

Comes with a silicone mat for safe placement of instruments. Case and silicone insert are autoclavable but should not be used for sterilizing surgical instruments.



### **Colored Plastic Instrument Cases**

This plastic case is compact, safe, secure, and slotted for placement of individual surgical instruments, aiding in the prevention of tip damage. Seven colors available. **Case is autoclavable but should not be used for sterilizing surgical instruments.** 



# SURGICAL ACCESSORIES

### **Roll Up Travel Case**

Vinyl dissection roll up case with Velcro strap to hold the case closed. Holds up to ten tools, one ruler holder, and one Velcro sealed pouch. Instruments not included.

36 x 23 cm (unfolded) Black Vinyl

No. 20831-10



### **Instrument Tip Protectors**

These tip protectors are a convenient, color-coded way to protect surgical instrument tips and edges. They provide exceptional protection and are vented to prevent rusting from condensation. Instrument tip protectors are recommended during handling, storage, and sterilization. Autoclavable.





1.6 x 19.1 mm

Pkg. of 25

Material: Plastisol

No. 29060-01

4.8 x 25.4 mm Pkg. of 25 Material: Plastisol

No. 29060-05





2 x 19 mm

Pkg. of 25

Material: Plastisol

No. 29060-02

Pka. of 25 Material: Plastisol

No. 29060-06









3.2 x 25.4 mm

Pkg. of 25

Material: Plastisol

No. 29060-04



Assorted pack, five of each size Pka. of 30 Material: Plastisol No. 29060-00

These vented tip protectors are designed to hold hemostats and rongeurs firmly while allowing the jaws, box hinge and tips to stay open for through sterilization. Autoclavable.

**Double Tip Protectors** 



12.7 x 2 mm Pkg. of 10 Material: Plastisol

No. 29066-01

9 x 25.5 mm Pkg. of 10 Material: Plastisol No. 29066-02



15 x 25 mm Pkg. of 10 Material: Plastisol

No. 29066-03



Pkg. of 10 Material: Plastisol

No. 29066-04



Assorted pack,

three of each size Pkg. of 12 Material: Plastisol

No. 29066-00

### **Burrs Stainless Steel**

Length: 44 mm Shaft Ø: 2.3 mm Pkg. of 10 burrs

$\mathop{Tip}_{(mm)} \emptyset$	Item
0.7 *	No. 19008-07
0.9 **	No. 19008-09
1.4	No. 19008-14

<sup>\*</sup> Use this size for screws No. 19010-10 (page 68). \*\* Use this size for screws No. 19010-00 (page 68).

### **Burrs Carbon Steel**

Length: 44 mm Shaft Ø: 2.3 mm Pkg. of 10 burrs



<sup>\*</sup> Use this size for screws No. 19010-10 (page 68).

### **Burrs Tungsten Carbide**

Length: 44 mm Shaft Ø: 2.3 mm Pkg. of 5 burrs



IIP Ø (mm)	Item
0.5	No. 19009-05
0.7 *	No. 19009-07
0.9 **	No. 19009-09

<sup>\*</sup> Use this size for screws No. 19010-10 (page 68). \*\* Use this size for screws No. 19010-00 (page 68).

### **Trephines Stainless Steel**

Side opening for clearing of material.

Length: 44 mm Shaft Ø: 2.3 mm Sold individually

Tip Outer Ø	Item
1.8	No. 18004-1
2.3	No. 18004-2
2.7	No. 18004-2
5	No. 18004-5



### **Circular Saw**

For cutting bone and other hard substances. Attaches to mandrel (No. 18000-18), on the right side.



Ø: 16 mm Thickness: 0.1 mm Pkg. of 10 saws

No. 18000-20

### Mandrel

Screw allows attachment of discs and circular saws for use with high speed drills.



Length: 45 mm Shaft Ø: 2.3 mm

No. 18000-18

F·S·T

<sup>\*\*</sup> Use this size for screws No. 19010-00 (page 68).

### **High Temperature Cautery Kit**

This high temperature cauterizer is a great alternative to disposable cautery units. Temperature at tips is approximately 1000 °C (1832 °F). Replacement cautery tips are available.







Plastic body Surgical Stainless Steel tips: 0.2 mm Ø Operation temperature: 1000 °C (1832 °F) Total weight: 70 g Length: 16 cm (without tip)

### **Cautery Kit**

### No. 18010-00

Cautery kit includes:

- 1 Cautery handle
- 2 AA batteries 2 Standard tips
- 2 Fine tips
- 1 Storage case



Replacement tips Single use

Straight / Standard Sterile: Yes Box of 10 tips

No. 18010-01



Replacement tips Single use Straight / Fine Sterile: Yes

Box of 10 tips

### No. 18010-02

### **Low Temperature Cautery Kit**

This low temperature cautery kit is a great alternative to disposable cautery units. Temperature at the tips is approximately 649 °C (1200 °F). Replacement cautery tips are available.





Plastic body Surgical Stainless Steel tips: 0.2 mm Ø Operation temperature: 649 °C (1200 °F) Total weight: 60 g Length: 12 cm (without tip)

### **Cautery Kit**

### No. 18019-00

Cautery kit includes:

- 1 Cautery handle
- 1 AA battery
- 2 Fine tips
- 2 Elongated fine tips
- 1 Storage case



Replacement tips Single use Straight / Fine Sterile: Yes

PRICES SUBJECT TO CHANGE PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER

Box of 10 tips No. 18019-01



Replacement tips Single use Elongated Straight / Fine Sterile: Yes

Box of 10 tips

No. 18019-02

### **Small Vessel Cauterizer Kit**

This small hand-held unit keeps bleeding to a minimum and therefore provides a surgical field clear of clamps and hemostats. The cauterizer should be activated momentarily so the delicate tips do not burn out. Temperature at tips is approximately 1000 °C (1832 °F). Replacement cautery tips are available. **Autoclavable.** 



Heavy gauge Stainless Steel body Platinum-iridium tips: 0.3 mm Ø Total weight: 100 g Length: 13.5 cm (without tip)

### **Cautery Kit** No. 18000-00

Cautery kit includes:

- 1 Cautery handle
- 2 AA Batteries
- 1 Straight fine tip
- 1 Angled fine tip
- 1 Storage case



Remove batteries before autoclaving.



Replacement tips Reusable Straight / Fine

No. 18000-01



Replacement tips Reusable Angled / Fine

No. 18000-02



Reusable Straight / Fine

No. 18000-03



Replacement tips Reusable Straight knob

No. 18000-04

# LABORATORY ACCESSORIES

### **Neuro Punches**

N.I.H. style neuro punches for punch removal of frozen brain microregions.



ID (mm)	OD (mm)	Gauge	Length (cm)	Item
0.31	0.61	24	4.2	No. 18036-24
0.41	0.7	22	4.5	No. 18036-22
0.69	1	19	4.5	No. 18036-19

### **Sample Corers**

For quantitative analysis. Comes with plunger.


8.5 cm



ID (mm)	OD (mm)	Gauge	Item
0.35	0.8	21	No. 18035-35
0.5	1	19	No. 18035-50
0.8	1.45	17	No. 18035-80
1	1.8	15	No. 18035-01
2	3	11	No. 18035-02
3	4	8	No. 18035-03
5	6	4	No. 18035-05

### **Absorbent Swabs Mounted**

These synthetic surgical sponges absorb fluids twice as fast as the best cotton sponges and have three times their fluid-holding capacity. They are made from strong, soft, lint-free biocompatible-proven material that is nonirritating to delicate tissue and leaves no fibrous residue that might cause wound contamination. The Absorption Spears are sterilized by radiation.



### **Micro Cleaning Brush**



### **Miniature File Set**

Comes with eight different file profiles in a vinyl pouch.





### **Sharpening Stones**

Carbon Steel

No. 30076-14

Natural Arkansas oil stones feature unequaled fineness, hardness and structural density. When used with mineral oil, they will produce a sharper cutting edge.

> 6 mm x 9 cm Square stone

> > No. 29008-01

6 mm x 9 cm Triangular stone

No. 29008-12

### 4-in-1 Screwdriver

With four hardened, interchangeable blades stored in handle: 1.5, 2, 2.5, and 3 mm.



Handle 9.5 cm, with attachment 11.5 cm Nickel handle No. 30051-10

### **Ceramic Blade Cutter**

Ceramic edge stays sharp 100 times longer than a steel edge. Rust and corrosion free. No trace metal contamination. Autoclavable.



Plastic handle with ceramic blade

Straight

13 cm Cutting Edge: 58 mm

No. 10025-70

Replacement blade

Straight

Cutting Edge: 58 mm

No. 10025-71

These scissors are extremely hard (54 on the Rockwell scale). Useful for cutting tin or aluminum sheets, bone, plastics, etc. The strong, 3 mm thick blades are magnetic. The plastic handles are anatomically shaped.





21 cm Stainless Steel (Inox) with plastic handle

No. 37500-00

### **Universal Scissors**

The scissor blades are made of high quality Stainless Steel and guarantee a perfect cut. Handles are made of glass fiber reinforced polyamide.






13 cm Stainless Steel (Inox) with plastic handle

No. 37500-05

### **Ruler with Gauge Measurements**

Front





Surgical Stainless Steel (Inox) 15.5 cm / 6 inch

No. 30085-15

### **Rulter with Metric Conversion Table**

Front



Back



50% of actual size

Surgical Stainless Steel (Inox) 30 cm / 12 inch

No. 30085-30

### **Ruler Label**

Used for stable measurements in all surgical platforms.






6 cm Pkg. of 25 labels

No. 29027-00

### **Plastic Vernier Caliper**

Tough polymer is unbreakable under normal conditions. Measures inside, outside, step and depth, metric and inch to 0.1 mm or 1/64". Maximum measurement: 153 mm.




**Digital Electronic Caliper** 

70% of

actual size

Wide screen LED display with the ability to measure in millimeters, inches and fractions. Measures inside, outside, step and depth. Battery, spare battery and plastic storage case included.



Castroviejo Surgical Caliper



Plastic 21 cm No. 30087-20 Surgical Stainless Steel (Inox) 23.5 cm

No. 30087-00

0-20 mm with 0.5 mm divisions Tip Ø: 0.5 mm

8.5 cm

No. 18000-35

### **Electronic Clock/Triple Timer**

Three activities can be timed up or down simultaneously or individually. Time can be displayed either in 12 or 24 hours and features a memory function. Maximum countdown time is 19 hours, 59 minutes, 59 seconds. Minimum time interval is 1 second. The timer housing includes a powerful magnet for mounting.

LCD display: 4.5 x 4.5 cm 7 x 8.5 x 1.5 cm AA battery (not included)

No. 27000-01



### **Gold Wire Contacts**

Matching pin and socket. Contacts are crimp terminated, but also can be used for solder connections.

Total length: 16 mm Ø:1.2 mm Opening for wire: 0.63 mm

Pkg. of 50 Male contacts

No. 19003-00

Pkg. of 50 Female contacts

No. 19003-01



F·S·T

### **Loupe with Micrometer and LED**

High quality loupe features a removable micrometer and ultra-bright LED lights surrounding the circumference for a defined image and clear contrast.



Magnification: 10x Field of view: 30 mm Working distance: 10 mm (without scale) Ultrabright LED

38 x 38 x 40 mm

No. 29026-10

### **Stage Micrometer**



76 x 26 mm 1 mm with 100 divisions (0.01 mm)

No. 29025-02

### **Headband Magnifiers**

Precision ground and polished optical glass lenses. The padded headband is lightweight and fully adjustable. Can be worn with spectacles.



Magnification: 2.5x Working distance: 20 cm (8 inch)

No. 28030-02

Magnification: 2.75x Working distance: 15 cm (6 inch)

No. 28030-03

Magnification: 3.5x Working distance: 10 cm (4 inch)

No. 28030-04

### Interchangeable Magnifier Spectsacles with LED

Economical and fully interchangeable headband magnifier with LED light.

MOST Popular

### Includes:

- Multiple Lens Magnifications (1.0x, 1.5x, 2.0x, 2.5x, 3.5x)
- LED Light Source
- Tilting Functionality
- Cushioned Nose Pads
- Adjustable Headband
- Interchangeable Lens Frames
- Requires 3 AAA Batteries (Batteries not included)



Magnification: 1x, 1.5x, 2x, 2.5x and 3.5x Field of view: Variable Working distance: Variable

No. 28080-10

### Galilean Flip-up

These loupes are lightweight and easy to use. Made from the finest optical coated lenses which deliver excellent clarity and high resolution. The loupes are mounted onto a strong yet flexible wraparound safety frame with a headband for stability and comfort. Ability to flip loupes up for use as safety glasses. Comes with mini screwdriver for fine adjustments, micro cloth and storage case.



Magnification: 3x Working distance: 40 cm (16 inch) Weight: 262 g

No. 28050-30

Magnification: 3.5x Working distance: 30 cm (12 inch) Weight: 77 g

No. 28050-35

### Eschenbach ridoMED

These magnifier spectacles allow full adjustments on almost all parts. Headband keeps spectacles steady and in place. Leather case and Headband included.



Magnification: 2.5x Working distance: 35 cm (14 inch) Weight: 97 g Plastic frame reinforced with fiberglass

No. 28040-30

### **Eschenbach maxDETAIL**

These spectacles are portable, lightweight, and easily adjustable. Designed with frosted objective lens housing which eliminates distracting reflections. Bendable ear pieces for a secure, comfortable fit. Protective case is included.



Magnification: 2x Field of view: 110 mm (4 inch) Working distance: 40 cm (16 inch) Weight: 114 g

No. 28030-20

maxDETAIL Spectacles shown with Eschenbach Clip-On LED Light (sold separately).

### **Eschenbach Clip-On LED Light**



No. 28030-25

### Eschenbach varioLEDflex

The varioLEDflex ensures ease and precision by its smooth flexible neck and tiltable magnifying lens. This magnifier's innovative features provide a spacious working platform suitable for any application. Mounting clamp provides uncompromised workspace due to its sturdy, slim design. The magnification lens is illuminated with shadow-free full or half ring optimum bright LEDs.

> Magnification: 2.5x Diopter: 6D Lens size: 76 mm Lens shape: Aspheric Working distance: Variable

Lens material: PXM® lightweight lens with patented Ceratec®scratch resistant coating

Illumination: SMD LED Power source: 120 V & 230 V

### No. 28015-25



### **Eschenbach Scribolux**

This magnifier provides an illuminated platform suitable for any working condition. The open design allows for a spacious working platform. Supplied in a microfiber pouch for protecting and cleaning the lens.

Magnification: 2.8x Diopter: 7/11D Lens size: 100 x 75 mm Lens shape: Aspheric Working distance: 7.5 cm (3 inch) Lens material: PXM® lightweight lens with patented Ceratec®scratch resistant coating Illumination: SMD LED Power source: 6 V (4 AA-Batteries included)

No. 28020-28



### **Lens Cleaning Kit**

This convenient kit contains everything you need to clean lens surfaces safely. Perfect for microscopes, camera lenses, binoculars, telescopes, graphic art optics and all other fragile lenses. Storage case included.

### Kit includes:

- 1 Bottle of lens cleaner (30 ml)
- 1 Ultra soft duster brush
- 1 Cleaning cloth (20 x 20 cm)
- 3 Books of lens paper (25 sheets each)



No. 29000-10

### **Lens Paper**

This premium grade optical lens tissue is extremely soft and safe for any glass or plastic optical lens surface.



15 x 10 cm Book of 50 sheets

No. 29015-50

Scalpel Handles  Forceps  207 – 212  Hemostats, Pin & Needle Holders, Retract Vessel Dilators & Clamps  213 – 217
Hemostats, Pin & Needle Holders, Retract Vessel Dilators & Clamps  213 – 217
Pin & Needle Holders, Retract Vessel Dilators & Clamps  213 – 217

### **FST® STUDENT LINE**



### **FEATURES**:

- Economy
- Quantity discounts
- Quality
- Less hand finishing
- Less expensive alloys

One of the hallmarks of FST has always been to accurately assess the needs of our customers. We gather information not only from the lab, but from conferences, journals and researchers. Whether it is designing new tools or finding new applications for existing instruments, we believe our best information comes directly from you.

Our customers wanted reduced costs and quantity discounts without sacrificing the precision and quality inherent in German instruments. To meet their needs, we have developed a full line of "student" instruments. They were created for use in the classroom where students are learning proper instrument care.

Unlike our competitors, the quality of these instruments does not drop dramatically with the price.

Our student line uses quality Surgical Stainless Steel and maintains the precision you expect from FST. In addition they carry our unconditional satisfaction guarantee while being offered at a significantly reduced price.

Not intended as a substitute for our highest quality instruments, the student line may require extra care as the alloys used are more susceptible to corrosion. Visit our website at **finescience.com** for more information.




Quantity discount for student instruments: 5-9=5%10+=10%Per item, selection cannot be assorted















11.5 cm

Cutting edge: 16 mm


No. 91461-11

Iris

MOST Popular



Sanvenero





for student 5-9 = **5%** 10 + = **10%** Per item, selection cannot be assorted for discount.

200

11.5 cm

Cutting edge: 16 mm

No. 91460-11

Wagner



5-9 = **5%** 10+= **10%** Per item, selection cannot be assorted





12 cm Cutting edge: 20 mm

No. 91468-12



**Spencer Ligature** 



9 cm Cutting edge: 16 mm No. 91476-09



Lister

Longer blunt blade prevents damaging tissue underneath skin without losing the cutting line.





17 cm

Cutting edge: 32 mm

No. 91410-17



Curved

Blunt/Blunt

17 cm

Cutting edge: 32 mm No. 91411-17

**Metzenbaum Baby** 

**Bone Scissors** 



Curved

Sharp/Sharp

Student Quality

10.5 cm

Cutting edge: 20 mm

No. 91405-10

**Strong Scissors** 



Quantity discount for student instruments:

5-9 = 5%

10 + = 10%

Per item, selection
cannot be assorted
for discount.

9 cm Cutting edge: 14 mm

No. 91604-09

Curved

Sharp/Sharp



Metzenbaum

Straight Blunt/Blunt 14.5 cm Cutting edge: 22 mm

No. 91416-14

Blunt/Blunt 18 cm Cutting edge: 24 mm

No. 91417-18

F·S·T **PRICES SUBJECT TO CHANGE** PLEASE SEE WEBSITE FOR CURRENT PRICING · ORDERING INFORMATION ON INSIDE COVER

Curved

Blunt/Blunt

14.5 cm

Cutting edge: 22 mm

No. 91419-14

for student

5-9 = 5% 10 + = 10% Per item, selection cannot be assorted for discount.

5-9 = **5%** 10+= **10%** Per item, selection cannot be assorted

for discount.

206

F·S·T

Scalpel Handle #3

3

For scalpel blades,

see page 78.

12 cm

No. 91003-12


Dumont #7

MOST POPULAR




Dumont #5









11 cm 0.17 x 0.1 mm Surgical Stainless Steel (Inox) No. 91197-00







Quantity discount for student instruments: 5-9 = 5% 10+ = 10% Per item, selection cannot be assorted for discount.



For scalpel blades,

see page 77.

13 cm


No. 91004-13

Straight

10.5 cm

Tip Ø: 0.3 mm

No. 91156-11



11 cm

0.2 x 0.2 mm

No. 91110-10

11 cm

0.2 x 0.1 mm

No. 91115-10

instruments: 5-9 = 5% 10 + = 10%

12 cm

0.2 x 0.2 mm

No. 91117-10

12 cm

0.2 x 0.2 mm

No. 91113-10

Straight

9 cm

Tip Ø: 0.75 mm

No. 91139-09

Straight


11.5 cm

Tip Ø: 1 mm

No. 91139-12

Micro Adson


Broad pattern with strong spring action.




Adson

Wide, strong grip with narrow tips.


Serrated Serrated 12 cm 12 cm Tip width: 0.7 mm Tip width: 1.5 mm No. 91106-12








Adson-Brown





**Adson with Cross** 

Serrations





Straight Serrated

Length (cm)	Item
11.5	No. 91100-12
13	No. 91100-13
14.5	No. 91100-14
16	No. 91100-16
17	No. 91100-17
	11.5 13 14.5 16

**Anatomical Narrow Pattern** 



Straight Serrated

Tip width	Length (cm)	Item
1.5	11.5	No. 91102-12
2.1	13	No. 91102-13
2.2	14.5	No. 91102-14
2.3	16	No. 91102-16





Straight 1x2 Teeth

Tip width	Length (cm)	Item
2.2	12	No. 91121-12
2.4	13	No. 91121-13
2.6	14.5	No. 91121-14
2.8	15.5	No. 91121-15
3	16	No. 91121-16
3.5	18	No. 91121-18
4	20	No. 91121-20
4.5	25	No. 91121-25

Surgical Narrow Pattern



Straight 1x2 Teeth

Tip width	Length (cm)	Item
2	11.5	No. 91123-12
2.2	14.5	No. 91123-14
2.4	16	No. 91123-16

No. 91119-12

5 - 9 = 5%

10 + = **10%** 

for discount.

Per item, selection cannot be assorted

No. 91127-12









Gillies



Serrated

15 cm

Tip dimension: 7.6 x 5.7 mm

No. 91126-15

Russian



Serrated

15 cm

Tip width: 1.2 mm


No. 91138-15


McIndoe





**Halsted-Mosquito** 







Hartmann



instruments:

5-9 = 5%

10 + = 10%

Per item, selection
cannot be assorted
for discount.

10 cm

No. 91303-10

Tip width: 1.4 mm

No. 91309-12

No. 91331-12

Curved Serrated

Tip width: 1.5 mm

Straight

14.5 cm

Tip width: 2.2 mm

No. 91108-14

10 + = **10%** Per item, selection cannot be assorted

for discount.

212

**Cross-Serrations** 

15.5 cm

Tip width: 2.1 mm

No. 91128-15

12.5 cm

Tip width: 1.4 mm

No. 91308-12

Crile

5-9 = **5%** 10+= **10%** 

for discount.

Per item, selection cannot be assorted



Kelly



**Rochester-Pean** 



Pean



DISCOUNT (%)

**Pin Holder** 

**Spring Loaded** 

No. 91706-05





Maximum spread: 1.1 cm

No. 91700-03

for student instruments: 5-9 = 5% 10+= 10%

### Dieffenbach



### Straight

35 mm Strong pressure

No. 91850-35



For suture needles, see

page 152 - 153.

Halsey

For pins, see page 145.

10 + = **10%** Per item, selection cannot be assorted for discount.

5 - 9 = 5%

### Nickel-plated

7.5 cm Jaw opening: 0 to 1 mm

No. 91606-07





Mayo-Hegar

Tip width: 2.6 mm

Clamping length: 21 mm Lockable: Yes Scissors: No

No. 91204-16



Senn-Miller



3 blunt teeth 7 x 8 mm and 1 blunt tooth 18 x 6 mm

No. 91791-16



**Vessel Dilator** 

11 cm Tip Ø: 0.3 mm

No. 91815-11



Per item, selection cannot be assorted for discount.

DISCOUNT (%)

00071-14	123	10035-00	76	10110-07	73	11003-13	109	11043-08	102	11200-33	91
00072-14	123	10035-05	76	10130-05	145	11003-14	109	11049-10	103	11203-23	91
00108-11	98	10035-10	76	10130-10	145	11003-16	109	11050-10	103	11203-25	91
00109-11	98	10035-12	76	10130-20	145	11003-18	109	11051-10	103	11210-10	91
00125-11	129	10035-15	76	10140-01	135	11003-20	109	11052-10	103	11210-20	91
00126-11	129	10035-20	76	10140-02	135	11006-12	107	11053-10	103	11215-01	61
00272-13	100	10035-25	76	10140-03	135	11006-15	101	11054-10	103	11215-02	61
00276-13	129	10035-30	76	10140-04	135	11008-13	110	11055-10	103	11220-21	90
00325-00	122	10050-00	75	10152-30	72	11008-15	110	11062-10	104	11223-20	90
00396-01	122	10052-11	75	10157-13	72	11009-13	110	11063-07	100	11231-20	88
00398-02	122	10053-09	75	10157-15	72	11009-15	110	11064-07	102	11231-30	88
00400-03	122	10053-13	75	10157-18	72	11010-17	115	11065-07	102	11241-30	90
00574-11	131	10055-12	70	10160-13	137	11011-17	115	11066-07	102	11242-40	90
00608-11	131	10056-12	70	10162-12	133	11012-18	115	11067-07	102	11251-10	85
00632-11	98	10058-15	71	10165-11	136	11012-21	115	11069-08	102	11251-20	85
00649-11	98	10059-13	70	10166-11	136	11012-25	115	11070-08	102	11251-21	86
10000-10	80	10059-15	71	10167-11	136	11013-18	115	11071-10	104	11251-23	85
10000-11	80	10060-00	80	10290-13	137	11013-21	115	11072-11	106	11251-30	85
10000-15	80	10060-13	80	10315-12	76	11014-18	114	11073-10	105	11251-33	87
10002-00	79	10061-12	135	10316-14	76	11015-18	114	11074-02	105	11251-35	87
10002-50	79	10062-12	135	10317-14	76	11016-17	115	11075-00	105	11252-00	83
10003-12	78	10063-15	134	10318-14	76	11017-17	115	11076-15	110	11252-20	85
10003-13	78	10064-14	135	10321-08	143	11018-12	107	11080-02	114	11252-21	86
10004-13	77	10065-15	134	10360-13	143	11019-12	107	11083-07	99	11252-23	85
10006-12	81	10066-15	134	10370-13	71	11020-12	107	11084-07	99	11252-30	85
10007-12	79	10073-14	70	10370-17	143	11021-12	112	11090-09	99	11252-40	85
10007-16	79	10075-16	69	10370-18	143	11021-14	112	11090-10	99	11252-50	83
10009-13	77	10080-05	69	10370-19	143	11021-15	112	11091-15	121	11253-10	86
10010-00	78	10081-10	69	11000-12	109	11022-12	107	11092-12	121	11253-20	84
10010-10	81	10082-15	69	11000-13	109	11022-14	112	11093-10	99	11253-21	84
10011-00	78	10083-20	69	11000-14	109	11022-15	112	11095-09	120	11253-25	84
10011-10	81	10087-12	139	11000-16	109	11023-10	112	11095-13	120	11253-27	87
10012-00	78	10088-15	137	11000-18	109	11023-12	112	11096-10	121	11253-29	84
10015-00	78	10089-11	139	11000-20	109	11023-14	112	11096-13	121	11254-20	83
10015-10	81	10090-13	138	11000-25	109	11023-15	112	11101-09	106	11255-20	85
10020-00	77	10090-17	138	11000-30	109	11024-12	110	11103-09	106	11258-20	92
10021-00	77	10091-12	139	11001-12	109	11024-18	115	11106-09	106	11260-20	88
10022-00	77	10092-12	140	11001-13	109	11025-18	115	11150-10	103	11270-20	89
10023-00	77	10093-13	140	11001-14	109	11026-15	113	11151-10	103	11271-30	89
10025-23	74	10094-13	140	11001-16	109	11027-12	107	11152-10	103	11272-30	89
10025-43	74	10095-12	69	11001-18	109	11028-15	113	11153-10	103	11272-40	85,89
10025-44	74	10096-16	69	11001-20	109	11029-14	112	11154-10	103	11272-50	89
10025-45	74	10098-06	141	11002-12	109	11030-12	108	11155-10	103	11273-20	84
10025-50	74	10099-12	142	11002-13	109	11031-15	113	11158-10	103	11273-22	84
10025-70	187	10099-15	142	11002-14	109	11032-12	107	11159-10	103	11274-20	89
10025-71	187	10100-30	74	11002-16	109	11036-20	115	11160-10	103	11280-11	131
10030-13	134	10100-45	74	11002-18	109	11036-90	115	11200-00	171	11281-11	131
10030-13	134	10110-01	73	11002-10	109	11040-08	102	11200-10	92	11282-11	131
10032-13	134	10110-03	73	11002-25	109	11041-08	102	11200-12	92	11293-00	88
10032-13	134	10110-05	73	11002-23	109	11041 08	102	11200-14	92	11294-00	90
.0000 10	107	10110 00	, ,	11000 12	100	11072 00	102	11200 17	52	11207 00	50

11295-00	85	11988-03	95	12500-12	149	14003-20	49	14062-10	37	14160-10	38
11295-10	85	11990-14	94	12501-13	149	14004-12	49	14063-09	37	14161-10	38
11295-20	83	11991-14	94	12502-12	150	14004-13	49	14064-11	44	14171-10	48
11295-51	85	11992-12	94	12502-14	150	14004-14	49	14065-11	44	14173-12	45
11297-00	89	11992-20	96	12503-15	151	14004-16	49	14066-11	44	14184-09	36
11297-10	89	12001-13	149	12504-15	151	14004-18	49	14067-11	44	14185-09	36
11370-31	97	12002-12	150	12510-14	148	14005-12	49	14068-12	39	14200-21	55
11370-32	97	12002-14	150	12512-12	150	14005-13	49	14069-12	39	14370-22	38
11370-40	97	12003-15	151	12565-14	148	14005-14	49	14070-12	39	14501-14	50
11370-42	97	12004-16	150	12875-13	146	14005-16	49	14071-12	39	14502-14	50
11373-12	97	12004-18	150	13002-10	120	14007-14	51	14072-10	47	14504-14	50
11373-22	97	12005-15	151	13003-10	120	14008-14	49	14074-09	43	14505-14	50
11390-25	98	12010-14	148	13004-14	118	14010-15	54	14074-11	44	14508-15	52
11393-04	100	12018-12	154	13005-14	118	14010-17	54	14075-09	43	14509-15	52
11393-05	100	12020-00	156	13006-12	116	14011-15	54	14075-11	44	14510-15	54
11399-80	97	12020-09	156	13007-12	116	14011-17	54	14076-09	46	14510-17	54
11399-87	97	12022-09	156	13008-12	117	14012-15	55	14077-09	44	14511-15	54
11412-11	95	12022-19	156	13009-12	117	14012-17	55	14077-10	44	14511-17	54
11413-11	95	12023-00	156	13010-12	117	14013-15	55	14078-10	42	14512-15	54
11445-12	95	12028-12	154	13011-12	117	14013-17	55	14079-10	42	14512-17	54
11480-11	96	12029-12	154	13012-12	117	14014-17	55	14080-11	47	14518-14	52
11485-11	96	12030-01	155	13013-14	119	14015-17	55	14081-09	46	14518-18	52
11487-11	96	12030-04	155	13014-14	119	14016-14	53	14082-09	46	14519-14	52
11506-12	107	12030-20	155	13015-14	119	14016-18	53	14083-08	35	14519-18	52
11515-02	61	12030-22	155	13018-14	118	14017-14	53	14084-08	35	14558-09	41
11602-14	109	12030-23	155	13019-14	118	14017-18	53	14084-09	36	14558-11	4′
11602-16	114	12035-00	155	13020-12	116	14018-13	52	14085-08	35	14559-09	4′
11608-15	110	12040-01	154	13021-12	116	14018-14	52	14085-09	36	14559-11	41
11616-15	108	12040-02	154	13108-12	117	14018-18	52	14086-09	47	14568-09	40
11617-12	108	12040-03	154	13109-12	117	14019-13	52	14088-10	38	14568-12	40
11626-11	87	12050-01	153	14000-12	49	14019-14	52	14090-09	36	14569-09	40
11627-12	108	12050-02	153	14000-13	49	14019-18	52	14090-11	36	14569-12	40
11650-10	103	12050-03	153	14000-14	49	14024-14	53	14091-09	36	14574-09	43
11651-10	103	12050-10	153	14000-16	49	14028-10	42	14091-11	36	14574-11	43
11652-10	103	12050-20	153	14000-18	49	14029-10	42	14094-11	37	14575-09	43
11700-00	105	12052-08	152	14000-20	49	14037-10	47	14095-11	37	14575-11	43
11700-01	105	12052-09	152	14001-12	49	14040-10	39	14101-14	51	14854-10	35
11700-02	105	12052-10	152	14001-13	49	14041-10	39	14105-12	45	14855-10	35
11700-03	105	12052-11	152	14001-14	49	14054-13	50	14106-09	48	14958-09	4′
11700-04	105	12060-01	146	14001-16	49	14058-09	40	14108-09	48	14958-11	4′
11851-20	93	12060-02	146	14001-18	49	14058-11	40	14110-15	54	14959-09	4′
11852-20	93	12061-01	146	14002-12	49	14059-09	40	14110-17	54	14959-11	4′
11911-11	93	12061-02	146	14002-13	49	14059-11	40	14111-15	54	15000-00	23
11915-11	93	12074-12	147	14002-14	49	14060-09	37	14111-17	54	15000-01	26
11917-12	93	12075-12	147	14002-16	49	14060-10	37	14116-14	53	15000-02	26
11923-13	94	12075-14	147	14003-12	49	14060-11	37	14116-18	53	15000-03	23
11972-20	96	12076-12	147	14003-13	49	14061-09	37	14117-14	53	15000-04	23
11980-13	95	12080-14	147	14003-14	49	14061-10	37	14117-18	53	15000-08	23
11988-01	95	12081-14	147	14003-16	49	14061-11	37	14130-17	50	15000-10	23
11988-02	95	12102-12	150	14003-18	49	14062-09	37	14131-14	51	15000-11	26

NUMERICAL INDEX

15000-12	26	15610-08	26	17004-05	159	18020-20	153	18058-15	125	18200-10	163
15001-08	23	15615-08	26	17005-03	158	18020-30	153	18058-18	125	18200-11	163
15002-08	23	15617-08	26	17005-04	159	18020-40	153	18059-14	124	18200-12	163
15003-08	24	15650-11	27	17006-05	161	18020-50	153	18060-18	165	18200-20	163
15004-08	24	15651-11	27	17007-08	161	18020-60	153	18060-20	165	18200-50	163
15005-08	24	15750-11	27	17008-07	160	18020-70	153	18060-22	165	18320-11	125
15006-09	28	15751-11	27	17009-07	160	18025-10	100	18060-24	165	18374-43	125
15007-08	25	15903-08	24	17009-08	160	18026-10	101	18061-10	165	18374-44	125
15008-08	24	15904-08	24	17010-07	160	18030-30	153	18061-15	166	18403-11	131
15009-08	25	16000-14	59	17011-10	162	18030-40	153	18061-18	166	18406-11	131
15010-09	29	16001-15	58	17012-11	162	18030-50	153	18061-20	166	18406-13	131
15010-10	29	16004-16	59	17012-13	162	18030-60	153	18061-22	165	18555-01	125
15011-12	31	16005-16	59	17020-17	157	18030-70	153	18061-24	165	18555-02	125
15012-12	31	16015-17	58	17021-13	157	18030-80	153	18061-38	165	18555-03	125
15013-12	31	16016-17	58	17022-13	157	18035-01	186	18061-50	165	18555-04	125
15015-11	30	16020-14	57	17023-13	157	18035-02	186	18061-75	165	18601-15	111
15017-10	28	16021-14	57	17300-02	159	18035-03	186	18062-12	128	18602-15	111
15018-10	30	16025-14	67	17300-03	159	18035-05	186	18062-13	128	18603-14	111
15019-10	30	16030-12	64	17300-04	159	18035-35	186	18063-13	128	18961-15	166
15020-15	32	16030-16	64	17370-53	159	18035-50	186	18063-20	165	18961-17	166
15021-15	32	16044-10	64	18000-00	185	18035-80	186	18064-20	165	18961-18	166
15023-10	27	16050-15	65	18000-01	185	18036-19	186	18065-20	165	18961-20	166
15024-10	27	16060-11	68	18000-02	185	18036-22	186	18067-10	165	18961-75	166
15025-10	27	16101-10	62	18000-03	185	18036-24	186	18067-11	133	18961-85	166
15030-14	34	16102-11	62	18000-04	185	18038-45	126	18067-20	165	19003-00	191
15031-14	34	16104-14	66	18000-10	130	18039-45	126	18067-50	165	19003-01	191
15033-09	28	16104-19	66	18000-18	183	18040-11	122	18067-75	165	19007-05	183
15036-14	34	16105-16	66	18000-20	183	18040-14	123	18080-01	127	19007-07	183
15037-14	34	16106-14	66	18000-30	169	18040-22	122	18080-02	127	19007-09	183
15040-11	29	16107-14	60	18000-35	190	18040-33	122	18080-03	127	19007-14	183
15042-11	29	16109-14	60	18000-45	174	18050-32	126	18090-45	175	19007-21	183
15043-08	25	16120-14	57	18000-46	174	18050-35	126	18090-46	175	19007-29	183
15044-08	25	16121-14	57	18000-47	175	18050-50	126	18090-50	175	19008-07	183
15070-08	25	16130-15	65	18000-50	174	18051-32	126	18090-51	175	19008-09	183
15071-08	25	16140-11	63	18000-51	174	18051-35	126	18105-01	186	19008-14	183
15110-08	33	16143-12	63	18000-52	175	18051-50	126	18105-04	186	19009-00	68
15111-08	33	16144-13	63	18004-18	183	18052-01	126	18131-12	129	19009-05	183
15123-12	32	16150-24	67	18004-23	183	18052-02	126	18132-12	129	19009-07	183
15124-12	32	16152-11	64	18004-27	183	18052-03	126	18144-30	125	19009-09	183
15200-00	33	16152-12	64	18004-50	183	18052-12	128	18155-13	104	19010-00	68
15370-50	29	16152-15	64	18010-00	184	18053-28	126	18196-10	186	19010-01	68
15370-52	29	16220-14	56	18010-01	184	18054-28	126	18200-01	163	19010-10	68
15371-92	30	16221-14	56	18010-02	184	18055-01	125	18200-02	163	19010-11	68
15396-00	26	17000-01	158	18019-00	184	18055-02	125	18200-03	163	20311-21	181
15396-01	26	17000-02	158	18019-01	184	18055-03	125	18200-04	163	20325-00	177
15400-12	32	17000-03	158	18019-02	184	18055-04	125	18200-05	163	20325-10	177
15401-12	32	17000-04	158	18020-00	153	18055-05	125	18200-06	163	20330-00	180
15403-08	24	17002-02	158	18020-01	153	18055-06	125	18200-07	163	20330-01	180
15404-08	24	17003-03	158	18020-02	153	18056-14	124	18200-08	163	20330-02	180
15514-12	31	17004-03	158	18020-03	153	18057-14	124	18200-09	163	20330-03	180

20330-04	180	26000-70	145	29027-00	189	91102-16	211	91400-12
20350-00	177	26001-25	145	29055-00	171	91106-12	210	91400-14
20810-01	178	26001-30	145	29060-00	182	91108-12	212	91401-10
20810-02	178	26001-35	145	29060-01	182	91108-14	212	91401-12
20830-00	181	26001-40	145	29060-02	182	91110-10	208	91401-13
20830-01	181	26001-45	145	29060-03	182	91113-10	208	91401-14
20830-02	181	26001-50	145	29060-04	182	91115-10	208	91402-10
20830-03	181	26001-55	145	29060-05	182	91117-10	208	91402-12
20830-04	181	26001-60	145	29060-06	182	91118-12	210	91402-13
20830-05	181	26001-65	145	29066-00	182	91119-12	210	91402-14
20830-06	181	26001-70	145	29066-01	182	91121-12	211	91404-12
20831-10	182	26002-10	145	29066-02	182	91121-13	211	91404-13
20840-01	179	26002-15	145	29066-03	182	91121-14	211	91404-14
20840-02	179	26002-20	145	29066-04	182	91121-15	211	91405-10
20840-03	179	26003-15	145	29070-01	172	91121-16	211	91405-11
20842-03	179	26007-01	145	29070-02	172	91121-18	211	91405-13
20850-00	178	26007-02	145	29070-03	172	91121-20	211	91405-14
20850-10	178	26007-03	145	29070-04	172	91121-25	211	91408-14
20850-11	178	26007-04	145	29070-05	172	91123-12	211	91409-15
24202-01	169	26007-05	145	29070-06	172	91123-14	211	91410-17
24202-02	169	26016-12	144	29070-07	172	91123-16	211	91411-17
24202-03	169	26018-17	144	29070-08	172	91126-15	212	91413-14
24202-04	169	26029-10	106	29070-09	172	91127-12	210	91413-18
24202-05	169	26030-10	106	29071-01	172	91128-15	212	91413-20
24202-06	169	27000-01	191	29071-02	172	91138-15	212	91416-14
24202-11	169	28000-06	173	29071-03	172	91139-09	209	91417-18
24202-12	169	28000-11	173	29071-04	172	91139-12	209	91419-14
24202-13	169	28000-14	173	29071-05	172	91150-20	207	91427-12
24202-14	169	28001-02	176	29071-06	172	91156-11	209	91460-11
24202-15	169	28015-15	176	29071-08	172	91162-12	210	91461-11
24202-16	169	28015-25	196	29071-09	172	91197-00	207	91468-12
24212-01	168	28020-20	176	29080-02	172	91201-13	216	91476-09
24212-02	168	28020-28	196	30051-10	187	91204-16	216	91500-09
24214-02	168	28030-02	192	30076-14	187	91303-10	213	91501-09
24215-02	167	28030-03	192	30085-15	189	91304-14	214	91604-09
24220-00	170	28030-04	192	30085-30	189	91304-16	214	91606-07
24221-10	170	28030-20	195	30087-00	190	91305-14	214	91700-03
24221-50	170	28030-25	195	30087-20	190	91305-16	214	91706-05
24221-90	170	28040-30	194	37500-00	188	91308-12	213	91790-16
24916-01	167	28050-30	194	37500-05	188	91309-12	213	91791-16
24916-02	167	28050-35	194	91003-12	206	91316-14	215	91815-11
26000-25	145	28080-10	193	91004-13	206	91316-16	215	91850-35
26000-30	145	29000-00	171	91100-12	211	91317-14	215	
26000-35	145	29000-10	197	91100-13	211	91317-16	215	
26000-40	145	29008-01	187	91100-14	211	91318-14	214	
26000-45	145	29008-12	187	91100-16	211	91319-14	214	
26000-50	145	29008-22	171	91100-17	211	91322-14	215	
26000-55	145	29015-50	197	91102-12	211	91323-14	215 213	
26000-60 26000-65	145 145	29025-02 29026-10	192 192	91102-13 91102-14	211 211	91331-12 91400-10	201	
20000-00	140	ZUZU-1U	134	J110Z-14	<b>411</b>	J1400-10	ZU I	

**ALPHABETICAL INDEX ALPHABETICAL INDEX** 

Absorption swabs	186
Adson DeBakey Forceps	108
Adson Forceps	107,210
Adson-Brown Forceps	108,210
Adson-Graefe Forceps	108
Agricola Retractor	159
Allis Clamp Forceps	121
Alm Retractors	160
Animal Identification	
Ear Punch	167-168
Ear Tags	170
Markers	169
Autoclave Tape	176

### В

D	
Baby-Mixter Hemostat	119
Backhaus Towel Clamp	120
Barraquer Retractors	217
BeeBee Bone Scissors	64
Biemer	
Clip Applicator	124
Clips	125
Bishop-Harman	102
Blair Retractor	157
Boehler Bone Cutter	65
Bone	
Austin Chisel	69
Dissecting Chisel	69
Forceps	68
Nipper	62
Scissors	205
Scraper	69
Screws	68
Spring Scissors	63
Wax	68
Bonn	
Forceps	99
Micro Probe	134
Scissors	35-36,47
Bonn-Strabismus Scissors	36
Breakable Scalpels	75
Bulldog Clamps	125

C	
Calipers	190
Cannulation System	130
Carbon Steel Blades	76
Castroviejo	
Needle Holder	146,148
Spring Scissors	27-28
Surgical Caliper	190
Cataract Knife	70
Cauterizer	184-185
Ceramic	
Scalpel Handle	74
Scalpels	74,187

Ceramic Coated	
Forceps	83,89
Iris Scissors	41
Spring Scissors	27
Circular Saw	183
Cleaning Brush	186
Cleveland Bone Cutter	66
Clips (Vascular)	
Biemer	125
Bulldog Artery Spring Clamp	125
Dieffenbach	217
Dietrich Bulldog	126
Double Clamps	122
Micro Serrefines	125
Mini Serrefines	126
S&T® Micro	122
Schwartz Micro Serrefines	126
Titanium	125
Clock (Electronic Timer)	191
Cohan-Vannas Spring Scissors	26
Colibri Retractors	158-159
Combination Spatula and Probe	138
Combination Spatula and Spoon	138
Cover Glass Forceps	105
Crile Hemostats	118,214
Crile-Wood Needle Holders	151
Cross-Action Forceps	105
Curettes	69
Cutters	61
Beebee Bone Scissors	64
Boehler Bone Cutter	65
Bone Nipper	62
Bone Pliers	67
Bone Scissors	64
Bone Shears	67
Bone Trimmers	60
Cartilage	60
Cleveland	66
Dumont	61
Liston Bone	66

Littauer Bone

Skull Flap

Malleus Nipper

D	
DeBakey Forceps	110
Decapitation Scissors	55
Diamond Coated Forceps	94
Diamond Knives	74
Dieffenbach Clamp	217
Dietrich Bulldog Clamps	126
Double Micro Clamp Applying Forceps	123
Double Micro Clamps	122
Dowell Spring Scissors	29
Drill Accessories	
Burrs	183
Circular Saw	183

Mandrel	183
Trephines	183
Dumont	61,83-92,131,207
Dumoxel	
Forceps	89

E	
Ear Punchs	167-168
Ear Tags	170
Episcleral Forceps	104
Eschenbach Magnifier	194-196
Clip-on LED light	195
Exclusive Line	
Forceps	93
Needle Holder	146
Rongeurs	56
Scissors	35

64

63

65

F	
Feeding Needles	
Reusable	165
Single Use	166
Feilchenfeld Forceps	209
File Set	187
Fomon Retractor	157
Forceps	95-96
Adson	107
Adson (Student)	210
Adson DeBakey	108
Adson-Brown	108
Adson-Brown (Student)	210
Adson-Graefe	108
Allis Clamp Forceps	121
Balanced	93,111
Biemer Clip Applicator	124
Bishop-Harman	102
Bone	68
Bonn	99
Bonn Micro	99
Ceramic Coated	83,89
Circle	106
Cover Glass	105
Cover Slip	87
Cross Action	105
DeBakey Round Handled	110
Diamond Coated	94
Double Micro Clamp Applying	123
Dumont	83-92,131
Dumont (Student)	207
Dumoxel	89_
Episcleral	104
Epoxy Coated	90-91
Feilchenfeld (Student)	209
Forceps Style Clip Applicator	124
Gerald	114
Gillies	113

Gillies (Student)

Insect Handling

McIndoe (Student) Micro Adson

Micro-Blunted Tips

Mirror Finish Moria Iris

Moria Paufique Moria Suture Tying

Moria Ultra Fine

Mouse Handling

Narrow Pattern

Organ Holding

Pierse Fixation

Replaceable Tips

Russian (Student)

S&T® SuperGrip

Self-Closing

Sliding Lock

Sorensen

Student

Taylor

Tissue

Titanium

Vessel Cannulation

Vessel Dilators (Student)

Forceps Style Clip Applicator

Forceps with Replaceable Tips

Friedmann-Pearson Rongeurs

Vessel Dilators

Friedmann Rongeurs

Suture Tying

Semken

S&T® Suture Tying

Semken (Student)

Splinter (Student) Standard Pattern

Standard Pattern (Student)

S&T® Micro Clamp Applying

S&T® Vessel Cannulation

Potts-Smith

Repair Kit

0-Rings

Octagon

Plastic

Ring

S&T®

Russian

Micro Adson (Student)

Laminectomy London

Gold Plated

Graefe

ing Needles	
Reusable	165
Single Use	166
nenfeld Forceps	209
Set	187
n Retractor	157
eps	95-96
Adson	107
Adson (Student)	210
Adson DeBakey	108
Adson-Brown	108
Adson-Brown (Student)	210
Adson-Graefe	108
Allis Clamp Forceps	121
Balanced	93,111
Biemer Clip Applicator	124
Bishop-Harman	102
Bone	68
Bonn	99
Bonn Micro	99
Ceramic Coated	83,89
Circle	106
Cover Glass	105
Cover Slip	87
Cross Action	105
DeBakey Round Handled	110
Diamond Coated	94
Oouble Micro Clamp Applying	123
Dumont	83-92,131
Dumont (Student)	207
Dumoxel	89
piscleral	104
poxy Coated	90-91
eilchenfeld (Student)	209
orceps Style Clip Applicator	124

### G

G	ialilean Loupes	19

Gavage Feeding Needles	165-166
GenoPunch	167
Gerald Forceps	114
Gillies Forceps	113,212
Gold Plated Forceps	93
Gold Wire Contacts	191
Goldstein Retractors	158
Graefe Forceps	102-103
Graefe Scissors	47
Grooved Director	137
Gross Anatomy Blade and Handle	80
Gross Anatomy Knives	72
Gross Anatomy Probe	137
Guthrie Retractor	157

### ш

212

93

90

114 212

107

210

84 85,94-95

97

98

100

97

115

109

171

102

113

104

105

115

171

95

106

113

212

98

123

98

100

131

96

110

212 87,104

111

209

109

211

208

115

112

131 129

217

101

124

95 59

56-57

89,114

99-101

102-103 106

Halsey Needle Holder	149,216
Halsted-Mosquito Hemostats	117,213
Hartmann Hemostats	120,213
Heavy Scissors	205
Heiss Retractor	162
Hemostats	116
Baby-Mixter	119
Crile	118
Crile (Student)	214
Halsted-Mosquito	117
Halsted-Mosquito (Student)	213
Hartmann	120
Hartmann (Student)	213
Kelly	118
Kelly (Student)	214
Mosquito	117
Pean (Student)	215
Providence-Hospital	119
Rochester-Pean (Student)	215
Ultra-fine	116
Hippocampal Dissecting Tool	142
Holders and Breakers	75
Hooks	134-135
Hot Bead Sterilizer	174-175

### Insect Handling Forceps

Insect Pins	14!
Insert Pins	14!
Instrument Care	
Double Tip Protectors	182
Forceps Repair Kit	17
Instrument Oil	17
Marking Tape	17:
Sharpening Stones	171,18
Tip Protectors	182
Instrument Cases	
Aluminum (Storage)	180-18
Plastic (Storage)	18
Sterilization	178-17

Instrument Oil	171
Iris Scissors	36-42,200
Iris Spatula	140

Kalt Needle Holders	147
Kaye Scissors	48
Kelly Hemostats	118,214

LaGrange Scissors	45
Large Loops Scissors	48
Lempert Rongeurs	59
Lens Paper	197
Lexer-Baby Scissors	42
Ligation Aids	128
Lister	
Scissors	51,203
Liston Bone Cutter	66
Littauer Bone Cutter	64
London Forceps	114
Lorna Towel Clamp	121
Luer Rongeurs	58

106

Magnifiers / Loupes	192-196
Cleaning Kit	197
Malleus Nipper	63
Mandrel	183
Markers	169
Marking Tape	172
Mathieu Needle Holder	148
Mayo Scissors	54,203
Mayo-Hegar Needle Holder	150,216
Mayo-Noble Scissors	55
Mayo-Stille Scissors	55
McIndoe Forceps	212
McPherson Forceps	100
Metzenbaum Baby Scissors	52,204
Metzenbaum Scissors	52-53,204
Michel Suture Clips & Applicators	154
Micro Clip Applicator	124
Micro Serrefines	125
Microsurgical Needles with Suture	Γhread
Attached	152
Mini Serrefines	126
Minutien Pins	145
Moria	
Bulldog Clamp	125
Forceps	97
Iris Scissors	38
Pascheff-Wolff Spring Scissors	30
Paufique Forceps	98
Retractors	159
Spoons	143

ALPHABETICAL INDEX

Spring Scissors	26
Suture Tying Forceps	100
Vannas-Wolff Spring Scissors	29
Vessel Clamp	125
Mosquito Hemostats Micro	117
Mouse Handling Forceps	115
Multipurpose Forceps	111

### N

Narrow Pattern Forceps	109
Needle Holders	
Castroviejo	146,148
Crile-Wood	151
Halsey	149
Halsey (Student)	216
Kalt	147
Mathieu	148
Mayo-Hegar	150
Mayo-Hegar (Student)	216
Olsen-Hegar	150
Round Handled	146
Tungsten Carbide	148-151
Needles	
Etched Tungsten Micro Needles	145
Holders	144
Holders (Student)	216
Insect	145
Insert	145
Minutien	145
Suture Needles	153
Neuro Punches	186
Noyes Spring Scissors	
Noyes opining ocissors	31
Nylon Suture	31 153

### 0

O-Rings for Forceps	171
Obsidian Stone Age Scalpels	73
Octagon Forceps	102
Olsen-Hegar Needle Holder	150
Optical Cleaning Kit	197
Organ Holding Forceps	113

### P

Г	
Pean Hemostats	215
Pearson Rongeurs	58
Pick	133
Pierse Fixation Forceps	104
Pin Holders	144,216
Plastic Forceps	105
Potts-Smith Forceps	115
Probes	135
Bonn Micro Probes	134
Double Pronged Pick	133
Providence-Hospital Hemostats	119
PTFE Feeding Needles	166

Repair Kit	171
Retraction System	163
Retractors	157-159,161,217
Agricola	159
Alm	160
Barraquer (Student)	217
Blair	157
Colibri	158-159
Fomon	157
Goldstein	158
Guthrie	157
Heiss	162
Senn-Miller (Student)	217
Weitlaner-Locktite	162
Reusable Feeding Needles	165
Ring Forceps	106
Rochester-Pean Hemostats	215
Rongeurs	
Friedmann	59
Friedmann-Pearson	56-57
Lempert	59
Luer	58
Pearson	58
Round Handled Needle Holde	rs 146-147
Rulers	189
Russian Forceps	113,212

S	
S&T®	
Clamp Applying Forceps	123
Forceps	98
Micro Clamps	122
SuperGrip Forceps	98
Suture Tying Forceps	100
Vessel Cannulation Forceps	131
Vessel Dilators	129
Sample Corers	186
Sanvenero Scissors	200
Scalpels	
Breakable	75
Cataract Knife	70
Ceramic Blade Cutter	187
Ceramic Blades	74
Circle Knife	71
Diamond Knives	74
Disposable Safety Scalpels	80
Holders and Breakers	75
Micro	76
Micro Knives	71,76
Safety Scalpel Blade Cartridges	81
Scalpel Blade Remover	79
Scalpel Blades	77-78
Scalpel Handles	77-79,81
Scalpel Handles (Student)	206
Schwartz Micro Serrefines	126
Scissors	
Bone	64

Bone (Student)	205
Bonn	35-36
Bonn Artery	47
Bonn-Gawenda	35
Bonn-Strabismus	36
Decapitation	55
Dissector	46
Graefe	47
Heavy (Student)	205
Iris	36-42
Iris (Student)	200
Kaye	48
LaGrange	45
Large Loops	48
Lexer-Baby	42
Ligature	46
Lister	51
Lister (Student)	203
Mayo	54
Mayo (Student)	203
Mayo-Noble	
Mayo-Stille	55
Metzenbaum	52-53
Metzenbaum (Student)	204
Metzenbaum Baby	52
Metzenbaum Baby (Student)	204
Micro	
Moria Iris	34
Sanvenero (Student) Shea	200
	45
Spencer (Ctudent)	46
Spencer (Student)	202 44
Stevens Tenotomy Strabismus	43
	43
Tissue Separating	
Universal Scissors	188
Utility Scissors	188
Wagner (Student)	39
Wagner (Student) Walton	202
	44
Screwdriver	187
Semken Forceps	110,212
Senn-Miller Retractors	217
Sharpening Stones	171,187
Shea Scissors	45
Sichel Knife	70
Silicone Mat	177
Silk Sutures	153
Single Use Feeding Needles	166
Skin Staplers & Applicators	155-156
Sliding Lock Forceps	104
Spatula	139,141
Combination Probe and Spatula	138
Combination Spoon and Spatula	138
Spencer Ligature Scissors	46,202
Spinal Cord Hook	133
Splinter Forceps	209
Spoons	143

### W

Wagner Scissors	39,202
Walton Scissors	44
Weitlaner-Locktite Retractor	162
Westcott Spring Scissors	30
Wraps	176

Spring Scissors	27,29,32
Bone Cutting	63
Castroviejo	27-28
Cohan-Vannas	26
Dowell	29
Micro	34
Moria	26,29-30
Noyes	31
Pascheff-Wolff	30
Titanium	26
ToughCut®	32
Trident II	33
Tungsten Carbide	31
Vannas	23-26,30,32
Vannas (Student)	200
Vannas-Tübingen	24
Vannas-Wolff	29
Westcott	30
Sterilization	
Containers	178-179
Disposable Wraps	176
Таре	176
Stevens Tenotomy Scissors	44
Storage Cases	180-181
Strabismus Scissors	43
Surgical Scissors	49-51
Student	201
Tungsten Carbide	40-41,43,50,52,54
Suture Needles	152-153
Suture Tying Forceps	99-101

### Т

Taylor Forceps	11
Timer (Electronic)	19
Tip Protectors	18
Tissue Forceps	11
Tissue Separating Scissors	4
Titanium	
Forceps	8
Titanium Forceps	11
Titanium Micro Point, Y-Hook, and Spatula	13
Towel Clamps	
Backhaus	12
Lorna	12
Travel Case for instruments	18
Trephines	18
Trident II Spring Scissors	3

### V

Vannas Spring Scissors	23-26	,30,32,200
Vannas-Tübingen Spring Sci	ssors	24
Vascular Occluders		127
Vessel Cannulation Forceps		131
Vessel Dilating Probe		137
Vessel Dilators		129,217
Vial Forceps		101
Von Graefe Knife		70

# **YOUR NOTES**

### FINE SCIENCE TOOLS AROUND THE WORLD

### SHIPPING GLOBALLY SINCE 1974



Most FST® instruments are made from Surgical Stainless Steel alloys with supreme German quality standards and premier craftsmanship. Most surgical instruments are shown in actual size. As these instruments are handcrafted, dimensions may vary slightly from those stated in this catalog.

The information in this catalog is accurate to the best of our knowledge. Designs, styles and dimensions depicted in this catalog are subject to change without notice. All recommendations or suggestions are made without guarantee, and we disclaim any liability incurred in connection with any recommendation. With the purchase of any product from this catalog the user assumes all risks and liabilities whatsoever in the use of such products. Neither seller nor manufacturer shall be liable in tort or contract for any loss or damage (direct, incidental, or consequential) arising out of the user, misses, or inability to use such product.

FST® is a registered trademark of Fine Science Tools (USA) Inc. and Fine Science Tools GmbH. ToughCut®, and CeramaCut® are registered trademarks of Fine Science Tools (USA), Inc.

© COPYRIGHT 2025/2026 The text, photos and drawings contained in this catalog are proprietary to Fine Science Tools Inc. and may only be copied, duplicated or otherwise reproduced with the written consent of Fine Science Tools Inc.



Scissors, Bone Instruments, Knives & Scalpels



### Fine Science Tools (USA) Inc.

4000 East 3<sup>rd</sup> Avenue, Suite 100 Foster City, CA 94404-4824

Phone +1 800 521 2109 / +1 650 349 1636 Fax +1 800 523 2109 / +1 650 349 3729 E-Mail america@finescience.com

Online finescience.com



Forceps/Tweezers, Hemostats & Vascular Instruments



### **Fine Science Tools GmbH**

Vangerowstraße 14 69115 Heidelberg

Germany **Phone** +49 62 21 - 90 50 50 Fax +49 62 21 - 90 50 590 E-Mail europe@finescience.com Online finescience.com



Probes, Needle Holders, Thread, Retractors & Clamps



### **Fine Science Tools**

Unit C, 10/F, Charmhill Centre 50 Hillwood Road, Tsimshatsui Hong Kong

Phone +852 6532 9229 E-Mail asia@finescience.com Online finescience.com



Surgical & Laboratory Accessories



Student Quality Instruments